

VWR Garments for Laboratories and Clean Rooms!

Savings!

Valid Until:
31 December 2018

Ref. Code:
VWR-GARM



Protection from Head to Toe - exclusively from Bio-Strategy!

A wide range of laboratory and clean room apparel.

Clean room garments are manufactured at special ISO accredited facilities. Data sheets and conformity declarations can be downloaded at uk.vwr.com. Stocked and on-demand garments available.


- Lab Coats
- Overall
- Aprons
- Isolation gowns
- Boot & Shoe covers
- Sleeve protectors
- Bouffant caps
- Beard covers
- Face Masks
- And much more!

Speak to your local Bio-Strategy Representative for further details or call Customer Service on 0800 34 24 66 or email sales.nz@bio-strategy.com



BIO-STRATEGY LIMITED

Tel: 0800 34 24 66
Email: sales.nz@bio-strategy.com
www.bio-strategy.com



Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.

DISPOSABLE LABORATORY APPAREL

Irradiated Coveralls, VWR

Irradiated coveralls are manufactured from a specially formulated breathable microporous fabric that provides significant fluid and barrier protection. Soft, cloth-like fabric. Coveralls are gamma irradiated to a Sterility Assurance Level (SAL) of 10⁻⁴ to address the concern of bio-burden contamination. Manufactured and packaged in an ISO Class 8 (FED-STD-209E Class 100,000/M3.5) cleanroom.

- Excellent fluid and particulate barrier,
- Superior breathability and water vapor transmission rate
- Extremely low level of particle shedding

Description	Pack Size	Availability	Code	List Price (excl. GST)	Sale Price (excl. GST)
Coveralls Irradiated with Attached Hood and Fluid-Impervious Boot Covers, Cleanroom ISO 8 M3.5					
Small	25	available	VWRU414004-449	\$319.00	\$234.00
Medium	25	available	VWRU414004-450	\$319.00	\$233.00
Large	25	available	VWRU414004-451	\$319.00	\$233.00
X-Large	25	available	VWRU414004-452	\$319.00	\$233.00
2X-Large	25	available	VWRU414004-453	\$319.00	\$233.00
3XL-Large	25	available	VWRU414004-454	\$319.00	\$233.00



ON SALE

Cleanroom overalls, sterile

These sterile single-use overalls are made of breathable SMS fabric. They are ideal for use in cleanroom production and research. Weight: 60 g, Individually double packed in easy-tear bags. 20 overalls per liner.

- With hood and two-way closure, elasticated wrist, ankles and back for tighter fit and easier donning, thumb loops (latex-free) and bound seams
- Suitable for use in ISO Class 5 + 6 cleanroom environments. Level I acc. IEST-RP-CC003.3
- Certifications: Type 5/6; Cat. III

Description	Pack Size	Availability	Code	List Price (excl. GST)	Sale Price (excl. GST)
Cleanroom overalls, sterile with hood, for ISO Class 5+6 cleanroom					
Small	20	enquire	VWRI113-7881	\$162.06	\$138.00
Medium	20	enquire	VWRI113-7882	\$169.83	\$144.00
Large	20	enquire	VWRI113-7883	\$163.17	\$139.00
X-Large	20	enquire	VWRI113-7884	\$180.93	\$154.00
2X-Large	20	enquire	VWRI113-7885	\$194.25	\$166.00
3XL-Large	20	enquire	VWRI113-7886	\$195.36	\$167.00



ON SALE



Cleanroom Lab coats, VWR Basic, SPP

Lightweight, and breathable spunbonded polypropylene (SPP) fabric, which protects against dirt, grime, and certain dry particulates in non-hazardous environments. Manufactured in an ISO certified facility.

- Tapered collar, inset sleeves, snap closure
- Elasticated or knitted cuffs, no pockets or three pockets
- Lightweight and breathable first line barrier protection



ON SALE

Description	Pack Size	Availability	Code	List Price (excl. GST)	Sale Price (excl. GST)
Cleanroom Lab coats, VWR Basic, SPP, disposable knitted cuffs, three pockets, white with collar					
Small	30	enquire	VWRI113-1111	\$91.67	\$78.00
Medium	30	enquire	VWRI113-1112	\$90.14	\$69.00
Large	30	enquire	VWRI113-1113	\$87.69	\$72.00
X-Large	30	enquire	VWRI113-1114	\$87.69	\$73.00
2X-Large	30	enquire	VWRI113-1115	\$95.46	\$85.50
3XL-Large	30	enquire	VWRI113-1116	\$96.57	\$80.50
Cleanroom Lab coats, VWR Basic, SPP, disposable knitted cuffs, elasticated cuffs, white, no pockets					
Small	25	enquire	VWRI113-1118	\$58.83	\$48.00
Medium	25	enquire	VWRI113-1119	\$63.10	\$48.50
Large	25	enquire	VWRI113-1120	\$65.04	\$51.50
X-Large	25	enquire	VWRI113-1121	\$62.39	\$54.50
2X-Large	25	enquire	VWRI113-1122	\$64.38	\$55.50

Lab coats, disposable, VWR Basic, SMS

SMS lab coats are manufactured from a latex-free, multilayer spunbonded-meltblown-spunbonded (SMS) fabric that offers high tensile strength while remaining soft, comfortable, and breathable.

- Tapered collar, inset sleeves, snap closure
- Knitted cuffs, three pockets
- Fluid and particulate barrier, light chemical splash resistance



ON SALE

Description	Pack Size	Availability	Code	List Price (excl. GST)	Sale Price (excl. GST)
Lab coats, disposable, VWR Basic, SMS, blue					
Small	30	enquire	VWRI113-1125	\$150.61	\$123.00
Medium	30	enquire	VWRI113-1126	\$138.04	\$123.00
Large	30	enquire	VWRI113-1127	\$144.95	\$127.00
X-Large	30	enquire	VWRI113-1128	\$144.34	\$127.00
2X-Large	30	enquire	VWRI113-1129	\$144.02	\$128.00
3XL-Large	30	enquire	VWRI113-1130	\$143.94	\$128.00

Laboratory coats, PP/PE, disposable

These disposable lab coats are made of PP (48g/m²) laminated with PE (24g/m²). They have two pockets (one breast pocket and one lower pocket).

- Length: 1020mm; 870mm (sleeves)
- Waterproof
- Tapered collar
- Bound seams with Press stud fastening
- Elasticated wrist



ON SALE

Description	Pack Size	Availability	Code	List Price (excl. GST)	Sale Price (excl. GST)
Laboratory coats, white PP/PE, waterproof, 2 pockets, disposable					
Large	30	enquire	VWRI113-8266	\$72.71	\$65.50
X-Large	30	enquire	VWRI113-8267	\$74.37	\$66.80
2X-Large	30	enquire	VWRI113-8268	\$74.37	\$66.80

Overalls, PP, with hood

Overalls made from non-woven PP. Ideal for non-critical applications.

- Hood, long sleeves and legs, Weight: 40g/m²
- Elasticated waist, wrist and ankles
- Bound seams and zip fastening

Description	Pack Size	Availability	Code	List Price (excl. GST)	Sale Price (excl. GST)
Overalls, PP, white, with hood					
Medium	25	enquire	VWRI113-8261	\$79.92	\$71.80
Large	25	enquire	VWRI113-8262	\$79.92	\$71.80
X-Large	25	enquire	VWRI113-8263	\$79.92	\$71.80




ON SALE

Aprons, LDPE/LLDPE

This disposable apron is made from a mixture of 65% LDPE and 35% LLDPE.

- Effective against dirt, moisture and oil
- Length: 1200 or 1400mm
- Colour: White

Description	Pack Size	Availability	Code	List Price (excl. GST)	Sale Price (excl. GST)
Aprons, disposable, LDPE/LLDPE					
700x1200mm	100	enquire	VWRI113-0740	\$16.76	\$14.20
<i>*replaced by VWRI113-8278 in November 2018</i>					



ON SALE

Isolation gowns SMS

These latex-free disposable gowns, with back fastening, are made of non-woven breathable SMS fabric. Ideal for use in healthcare environments and hospitals.

- Long sleeves, knitted cuffs, Round neck
- Ties on neck and waist
- Weight: 45 or 60g/m²

Description	Pack Size	Availability	Code	List Price (excl. GST)	Sale Price (excl. GST)
Isolation gowns SMS, disposable, 115x136 cm, back fastening, VWR					
45g/m ² White	50	enquire	VWRI113-8315	\$118.66	\$99.50
45g/m ² Blue	50	enquire	VWRI113-8316	\$118.66	\$109.00
60g/m ² White	50	enquire	VWRI113-8321	\$145.82	\$125.00
60g/m ² Blue	50	enquire	VWRI113-8322	\$142.37	\$132.00



ON SALE

Boot covers, non-skid, VWR Maximum, SF

VWR Maximum SF boot covers are the cleanest and most durable within the entire VWR clothing line. Made from spunbonded polypropylene fabric (SPP) coated with breathable polyethylene film, these garments have been designed to optimise barrier protection, cleanliness, and comfort properties. Length: 356mm (universal size, other sizes on demand)

- Maximum skid resistance on wet and smooth surfaces
- Clean, low linting fabric, Fluid impervious material
- Tested and manufactured in an ISO certified facility
- Lot controlled and validated through independent laboratory testing

Description	Pack Size	Availability	Code	List Price (excl. GST)	Sale Price (excl. GST)
Boot covers, non-skid, VWR Maximum, SF, white, with ankle ties, clean room ISO certified					
Universal size	100 pairs	enquire	VWRI113-1190	\$349.00	\$269.00



ON SALE

Shoe covers, non-skid, VWR Maximum SF

VWR Maximum SF shoe covers are the cleanest and most durable within the entire VWR clothing line. Made from spunbonded polypropylene fabric (SPP) coated with breathable polyethylene film, these garments have been designed to optimise barrier protection, cleanliness, and comfort properties. Length: 359mm (universal size, other sizes on demand)

- Maximum skid resistance on wet and smooth surfaces
- Clean, low linting fabric, Fluid impervious material
- Tested and manufactured in an ISO certified facility
- Lot controlled and validated through independent laboratory testing

Description	Pack Size	Availability	Code	List Price (excl. GST)	Sale Price (excl. GST)
Shoe covers, non-skid, VWR Maximum SF, white, ISO certified					
Universal size 359mm	200	enquire	VWRI113-1369	\$268.62	\$206.00



Shoe covers, CPE

Basic shoe covers manufactured from a low-density chlorinated polyethylene fabric, which is designed to protect against dirt, grime, and liquid splash.

- CPE is water resistant and provides good resistance against chemicals
- Economical and lint-free, one size fits all

Description	Pack Size	Availability	Code	List Price (excl. GST)	Sale Price (excl. GST)
Shoe covers, CPE					
Universal size, white	1000	enquire	VWRI113-1491	\$158.73	\$94.50
Universal size, blue	1000	available	VWRI113-1492	\$158.73	\$94.50



Cleanroom sleeve protectors VWR Advanced, SF

VWR Advanced SF sleeve protectors are manufactured from spunbonded polypropylene fabric (SPP) coated with breathable polyethylene film.

- Elasticated cuffs, with sonically welded seams, soft cloth-like fabric
- Durable fabric with good breathability and water transmission rate, Significant fluid and particle barrier
- Particle shedding: Level II Helmke drum classification
- Tested and manufactured in an ISO certified facility under stringent process controls
- Lot controlled and validated through independent laboratory testing

Description	Pack Size	Availability	Code	List Price (excl. GST)	Sale Price (excl. GST)
Cleanroom sleeve protectors, elasticated cuffs, VWR Advanced SF, ISO certified					
495mm	300	enquire	VWRI113-1194	\$137.00	\$105.00



Bouffant caps, VWR® Basic

These disposable caps are made from a non-woven, spunbonded polypropylene (SPP) fabric, which protects against dirt, grime and certain dry particulates in nonhazardous environments. Weight: 10g/m²

- Free from natural rubber latex
- Available in a wide range of colours

Description	Pack Size	Availability	Code	List Price (excl. GST)	Sale Price (excl. GST)
Bouffant caps SPP, latex free, white, VWR Basic					
530mm	1000	enquire	VWRI113-1459	\$53.10	\$51.80



Beard covers, VWR Basic, SPP

The VWR Basic SPP beard covers are manufactured from a strong, lightweight, and breathable spunbonded polypropylene (SPP) fabric, which protects against dirt, grime, and certain dry particulates in nonhazardous environments. Diameter (stretched): 406mm

- Lightweight and breathable first line barrier protection
- Latex-free version

This garment is ideal for less critical areas or pre-gowning entry rooms. It is tested and manufactured in an ISO certified facility, lot controlled and validated through independent laboratory testing.

Description	Pack Size	Availability	Code	List Price (excl. GST)	Sale Price (excl. GST)
Beard covers SPP, latex free, white, ISO certified, VWR Basic					
406mm	1000	available	VWRI113-1464	\$67.27	\$59.50



ON SALE

Hygiene masks, 2-ply non woven PP, VWR®

This hygiene mask, with nose bridge, is made from 2 ply, non woven PP.

- Density: 20g/m² (front layer); 25g/m² (back layer)
- Latex-free
- Available with ear loops or ties (optional)

Description	Pack Size	Availability	Code	List Price (excl. GST)	Sale Price (excl. GST)
Hygiene masks, 2-ply non woven PP, with earloops, VWR					
White	50	enquire	VWRI113-0713	\$5.80	\$4.10



ON SALE

Hygiene masks, 3 ply non-woven PP, VWR®

Hygiene masks made from 3 ply, non woven PP featuring a nose bridge and ear loops or ties.

- Particle size: 0.1 µm and density: 20g/m² (front layer); 25g/m² (back layer); 20g/m² (centre)
- Latex-free
- Bacterial filtration efficiency: 90% and particle filtration efficiency: 98.5%

Description	Pack Size	Availability	Code	List Price (excl. GST)	Sale Price (excl. GST)
Hygiene masks, 3 ply non woven PP with earloops, VWR					
Green	50	enquire	VWRI113-0714	\$7.20	\$5.10
Blue	50	enquire	VWRI113-0716	\$5.84	\$5.10



ON SALE

Cleanroom face masks, Maximum

Maximum cleanroom masks feature three ply construction for superior particle and bacterial filtration efficiency. These masks are latex-free. 3-ply masks provide high bacterial and particle filtration efficiency, 50 masks per cleanroom compatible, heat-sealed bag; 10 bags per case.

- Highly breathable and fluid resistant, moisture resistant, sterile mask features a unique, lightweight inner layer of film for enhanced moisture resistance
- A full length plastic arch holds the fabric away from the wearer's mouth to allow cooler, more comfortable breathing, malleable noseband ensures a secure seal
- Particle shedding (particles: ≥0,5 µm/cu. ft.): 196 (IEST-RP-CC003.3), bacterial filtration efficiency: >99 % (ASTM F2101), particle filtration efficiency: > 99 % (ASTM F2299), Delta P breathability (mm H₂O/cm²): 1.9 - 2.2 (MIL-M-36954C)
- Rigorously tested and manufactured in an ISO certified facility; products are validated through independent laboratory testing

Description	Pack Size	Availability	Code	List Price (excl. GST)	Sale Price (excl. GST)
Cleanroom face masks, Maximum, 3-ply, ISO certified, 175x95mm, VWR					
With earloops, blue	500	enquire	VWRI113-6220	\$75.44	\$69.50
With ties, white	500	enquire	VWRI113-6221	\$71.00	\$69.00



ON SALE

Laboratory coats, ladies and men, white cotton, VWR

These reusable coats for women are either made of 100% cotton or 65% polyester/35% cotton (optional). Both fabrics are robust, tear resistant and durable. The coats are available with a tailored or stand-up collar. Density: 215g/m², Colour: White

- Concealed press stud fastening and two-step adjustable cuffs
- Two side pockets and one breast pocket
- Fitted style, with back belt and side slit
- Washable at 95°C, suitable for industrial laundering

Description	Pack Size	Availability	Code	List Price (excl. GST)	Sale Price (excl. GST)
Laboratory coat, ladies, white cotton, tailored collar					
Size 8	1	enquire	VWRI113-1682	\$56.25	\$41.50
Size 10	1	enquire	VWRI113-1683	\$56.25	\$41.50
Size 12	1	enquire	VWRI113-1684	\$43.32	\$41.50
Size 14	1	available	VWRI113-1685	\$43.32	\$41.50
Size 16	1	enquire	VWRI113-1686	\$56.31	\$41.50
Size 18	1	enquire	VWRI113-1687	\$56.25	\$41.50
Size 20	1	enquire	VWRI113-1688	\$56.25	\$41.50
Size 22	1	enquire	VWRI113-1689	\$56.25	\$41.50
Size 24	1	enquire	VWRI113-1690	\$56.25	\$41.50
Size 26	1	enquire	VWRI113-1691	\$56.25	\$41.50
Laboratory coat, men, white cotton, tailored collar					
Size 36	1	enquire	VWRI113-1725	\$43.32	\$41.50
Size 38	1	enquire	VWRI113-1726	\$56.25	\$41.50
Size 40	1	enquire	VWRI113-1727	\$56.25	\$41.50
Size 42	1	enquire	VWRI113-1728	\$56.25	\$41.50
Size 44	1	enquire	VWRI113-1729	\$56.25	\$41.50
Size 46	1	enquire	VWRI113-1730	\$56.25	\$41.50

Also available: men, ladies coats with stand up collar and additional sizes, enquire today!



**ON
SALE**



Bio-Strategy Ltd.

T: 0800 34 24 66
E: sales.nz@bio-strategy.com
www.bio-strategy.com



bio-strategy
delivering technology