

Huge savings on VWR Next Generation pipette tips!



Valid Until:
31 December 2018

Ref. Code:
VWRI-NEXT

SAVE
NOW



Request a copy of the compatibility guide with further product details and the full range!

Ask us for samples or a demonstration and start saving!

Save money and reduce waste with VWR Next Generation Tips!

VWR Pipette tips NEXT GENERATION refill packs, racks and refill boxes

These tips are made of high quality 100% pure virgin medical grade PP. Made in the highest quality moulds, these tips can be used whenever the work demands a higher standard.

- A wider seal ensures that tips will seal easily and consistently without leaking
- Guaranteed free of detectable DNA, DNase/RNase, endotoxins and heavy metals
- Sterile products are tested for endotoxin (pyrogen) contamination and are certified to USP standards, Resins are pre-tested for metal contamination
- Low carbon footprint makes them an ideal solution to waste reduction goals for the laboratory
- Refill racks without having to touch the new tips

Tip styles:

Bevel Point™: Designed to retain less fluid and are a good choice for a multitude of procedures

UltraFine™: Thin flexible section at the end that is non-bevelled and very flexible. This helps you dispense the last tiny amount of sample from the tip.

FlexTop™: Uniquely thin and flexible top that conforms to the pipette barrel unlike any other tips. Low force ejection help reduce the potential for repetitive strain injuries.

Low Retention: Using special resin that reduces the potential for sample retention inside the tip. Ideal for protein research. Lot-certified and protease/nuclease free.

Sterile Filter tips, SLIK, Low Retention

Low retention tips offering unsurpassed fluid retention performance and proven aerosol blocking with the endotoxin-free, 10-micron porous hydrophobic filter. Ideal for PCR and DNA manipulation.

- Transparent PP resins, with PE filter, sterile
- Up to 10 times less fluid retention than standard tips
- Sample can be completely recovered without cutting into tip
- Lot-certified purity (free from RNase/DNase, human DNA and endotoxin, non-pyrogenic)
- Tips and racks are autoclavable to 121°C for maximum 1 min (dry cycle is not recommended)

Rainin® LTS® style Pipette Tips

The economic alternative for your LTS® pipettor in bulk or Next Generation Stacks

Speak to your local Bio-Strategy Representative for further details or call Customer Service on 0800 34 24 66 or email sales.nz@bio-strategy.com

View codes and pricing on the next page...

BIO-STRATEGY LIMITED

Tel: 0800 34 24 66
Email: sales.nz@bio-strategy.com
www.bio-strategy.com



Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.

**SAVE
NOW**

Description	Length	Volume	Pack Size	Availability	Code	Sale Price (excl. GST)
Next Generation Refill System						
1250µl UltraFine™, FlexTop™, extended, Next Generation refill (5 refills, 96 each)	88.90mm	100-1250µl	480	available	VWRI613-0738	\$22.20
300µl UltraFine™, FlexTop™, Next Generation (6 refills, 96 each)	71.12mm	1-300µl	576	on demand	VWRI613-0723	\$23.80
200µl Bevel Point™ (yellow) Next Generation refill (10 refills, 96 each)	49.66mm	1-200µl	960	available	VWRI613-0732	\$17.90
10µl Extended length, Tubegard™ ring tips, Next Generation refill (14 refills, 96 each)	38.10mm	0.1-10µl	1344	available	VWRI613-0735	\$31.15
10µl UltraFine™, Next Generation refill, (13 refills, 96 each)	31.42mm	0.1-10µl	1248	available	VWRI613-1646	\$28.63
Next Generation Refill Rack Boxes (empty)						
Refillable boxes for VWR Next Generation refills for 1-250µl tips, empty			10	available	VWRI613-0307	\$11.20
Refillable boxes for VWR Next Generation refills for 100-1250µl tips, empty			6	available	VWRI613-0306	\$11.40
Bulk Tips						
5ml Standard tips, bulk	144.27mm	1000-5000µl	250	available	VWRI613-0830	\$25.10
5ml Standard tips (for Eppendorf/Sartorius), bulk	123.06mm	1000-5000µl	250	available	VWRI613-0339	\$19.70
10ml Standard tips, bulk	155.70mm	5000-10000µl	250	on demand	VWRI613-0831	\$64.00
Filter Tips						
10µl Filter tips, low retention SLIK, Sterile, 12 racks, 96 each, VWR	31.42mm	10µl	1152	available	VWRI732-1486	\$61.30
10µl Extended Filter tips, low retention SLIK, Sterile, 10 racks, 96 each, VWR	38.10mm	10µl	960	available	VWRI732-1487	\$61.30
10µl Extra long, Low retention tips, sterile, 10 racks, 96 each, VWR		10µl	960	available	VWRI732-1614	\$61.30
20µl Filter tips, low retention SLIK, Sterile, 10 racks, 96 each, VWR	49.66mm	20µl	960	available	VWRI732-1488	\$61.30
200µl Filter tips, low retention SLIK, Sterile, 10 racks, 96 each, VWR	52.07mm	200µl	960	available	VWRI732-1489	\$61.30
1000µl Filter tips, low retention SLIK, Sterile, 6 racks, 96 each, VWR	88.90mm	1000µl	576	available	VWRI732-1491	\$62.40
Rainin® LTS® style Pipette Tips						
20µl Rainin® LTS® tips, bulk	48.26mm	0.5-20µl	1000	available	VWRI613-0578	\$18.29
200µl Rainin® LTS® tips, bulk	51.31mm	1-200µl	1000	available	VWRI613-0581	\$21.41
1000µl Rainin® LTS® tips, bulk	88.90mm	100-1000µl	1000	available	VWRI613-0584	\$24.23
20µl Pipette tips for Rainin® LTS® pipettes, Next Generation refill (11 refills, 96 each), sterile	48.26mm	0.2-20µl	1056	on demand	VWRI613-0722	\$42.58
250µl Rainin® LTS® tips, Next Generation refill (10 refills, 96 each)	51.31mm	0.5-250µl	960	on demand	VWRI613-0727	\$36.00
1000µl Rainin® LTS® tips, Next Generation refill (5 refills, 96 each)	88.90mm	100-1000µl	480	available	VWRI613-0729	\$26.29



BIO-STRATEGY LIMITED

Tel: 0800 34 24 66
Email: sales.nz@bio-strategy.com
www.bio-strategy.com