

## Axygen® Long Reach Pipette Tips

**Special pricing - 30% discount, samples available on request**

### Ideal for pipetting from long and narrow tubes

Axygen's NEW Long Reach Tips are longer and thinner to help you minimise contamination when aspirating from long and narrow storage tubes, cell culture flasks, plates and sample bags. The tips are elongated, clear and compatible with the vast majority of research grade pipetters.

- Slender Elongated Tip
- Graduation marks for easy visualisation of liquid volume
- DNase, RNase, pyrogen free
- Autoclavable (121°C for 15 minutes) or purchase already sterilised tips
- 10µl tip length up to 45mm
- 200µl tip length up to 77.5mm
- 1000µl tip extended length (102mm)


**Speak to your local Bio-Strategy Representative for further details or call Customer Service on 0800 34 24 or email [sales.nz@bio-strategy.com](mailto:sales.nz@bio-strategy.com)**

**VALID UNTIL:  
31 March 2019**

**Ref. Code:  
AXYG-LRTIPS**

**30%  
OFF**

Description	Pack Size	Availability	Code	List Price (excl. GST)	Sale Price (excl. GST)
Axygen 10µl long reach tips, bulk, non-sterile	1000 items	stocked	AXYGT-10XT-C	<del>\$49.15</del>	<b>\$34.41</b>
Axygen 10µl long reach tips, rack, sterile	960 items	stocked	AXYGT-10XT-C-R-S	<del>\$72.54</del>	<b>\$50.78</b>
Axygen 200µl long reach tips, bulk, non-sterile	1000 items	stocked	AXYGT-200XT-C	<del>\$52.54</del>	<b>\$36.78</b>
Axygen 200µl long reach tips, rack, sterile	960 items	stocked	AXYGT-200XT-C-R-S	<del>\$88.14</del>	<b>\$61.70</b>
Axygen 1000µl long reach tips, bulk, non-sterile	768 items	stocked	AXYGT-1000XT-C	<del>\$48.13</del>	<b>\$33.69</b>
Axygen 1000µl long reach tips, rack, sterile	768 items	stocked	AXYGT-1000XT-C-R-S	<del>\$79.33</del>	<b>\$55.53</b>



Bio-Strategy Ltd.

T: 0800 34 24 66

E: [sales.nz@bio-strategy.com](mailto:sales.nz@bio-strategy.com)

[www.bio-strategy.com](http://www.bio-strategy.com)

<https://shop.bio-strategy.com>

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.

# Axygen® Long Reach Pipet Tips

Ideal for pipetting from long and narrow tubes

**AXYGEN**®

A Corning Brand


**NEW**

When pipetting, are you currently disturbed because you are unable to aspirate precious reagents/samples from long and narrow storage tubes? Does it bother you that you might contaminate the micropipettor during pipetting?

Your help is here. Axygen® Long Reach Tips are elongated, clear, and compatible with vast majority of research grade pipettor brands. The tips are suitable for basic research experiments, and also able to meet the needs of high-end precision pipetting.

## Features

- ▶ Slender Elongated Tip
- ▶ Graduation marks for easy visualization of liquid volume
- ▶ DNase, RNase, pyrogen free
- ▶ Autoclavable (121 °C, 15 minutes)
- ▶ Available in various packaging sizes
- ▶ Products are manufactured in a facility that conforms to ISO14644-1 class ISO 8 qualification
- ▶ Materials conforming with US FDA USP Class VI level
- ▶ Resin is delivered by central control feeding system during the manufacturing process, hence eliminating external contamination



Using the 10µl long Reach Tips, you can pipette antibody (100µl / tube, 1mL tube) directly from the bottom of the tube. Additionally, because the tips are longer and thinner with a smaller surface area, this reduces the volume of residual antibody that sticks on the outside of the pipette tip compared to using normal tips.

The length of 200µl long reach tip ensures that the micropipettor shaft and tip ejector arm is not contaminated with dyes during buffer loading or cell staining procedures.

If you work with Matrigel, we recommend the use of 200µl or 1mL extended tips. The tips can reach the bottom of the Matrigel bottle. The shaft and ejector arm of the pipettor will not come into contact with the mouth of the bottle and hence reducing the chances of contamination.




Figure 1. Changing the media of T25 cell culture flask with a 1 mL Long Reach Tip is simple, efficient and convenient



Figure 2. Micro-pipetting with 10µl Long Reach Tips for precision pipetting


## 10 µl Long Reach Pipet Tips

- Up to 45mm, longer than standard pipette tips, suitable for micro-pipetting, allows easy access to the bottom of tubes
- Graduation marks at 2µl, 5µl, and 10µl
- Volume Range: 0.1µl-10µl
- Universal fit and compatible with majority of micropipettor brands




## 200 µl Long Reach Pipet Tips Coming soon

- Up to 77.5mm, you can easily reach the bottom of the container, without damaging the pipetting pipette. 10µl, 50µl, 100µl and 200µl Minimum fluid volume 1 µl Compatible with the vast majority of pipettes
- Graduation marks at 2µl, 5µl, and 10µl
- Volume Range: 1µl-200µl
- Universal fit and compatible with majority of micropipettor brands



## 1000 µl Long Reach Pipet Tips

- Extended length (102mm) reaches the bottom of tubes without contaminating the pipettor shaft
- 200µl, 500µl and 1000µl mark at the mark Compatible with the vast majority of pipettes Deep hole containers are most suitable
- Universal fit and compatible with majority of micropipettor brands
- Ideal for pipetting from deep vessels



## Ordering Information

CAT.NO.	Description	Packaging Format
T-10XT-C	10µl long reach tips, Clear, Bulk Pack, Nonsterile	1000 Tips/Pack, 10 Packs/Case
T-10XT-C-R-S	10µl long reach tips, Clear, Rack Pack, Sterile	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
T-200XT-C	200µl long reach tips, Clear, Bulk Pack, Nonsterile	1000 Tips/Pack, 10 Packs/Case
T-200XT-C-R-S	200µl long reach tips, Clear, Rack Pack, Sterile	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
T-1000XT-C	1000µl long reach tips, Clear, Bulk Pack, Nonsterile	768 Tips/Pack, 10 Packs/Case
T-1000XT-C-R-S	1000µl long reach tips, Clear, Rack Pack, Sterile	96 Tips/Rack, 8 Racks/Pack, 5 Packs/Case

### Corning Incorporated Life Sciences

836 North St.  
Building 300, Suite 3401  
Tewksbury, MA 01876  
t 800.492.1110  
t 978.442.2200  
f 978.442.2476

[www.corning.com/lifesciences](http://www.corning.com/lifesciences)

### Technical Support Hotline

Singapore: +65 6622 3609  
Australia: +61 2 8098 9173  
New Zealand: +64 9913 5888

[ScientificSupportAP@corning.com](mailto:ScientificSupportAP@corning.com)

**ASIA/PACIFIC**  
Australia/New Zealand  
t 61 427286832

**Korea**  
t 82 2-796-9500  
f 82 2-796-9300

**Singapore**  
t 65 6572-9740  
f 65 6861-2913

**India**  
t 91 124 4604000  
f 91 124 4604099

**CORNING**

**FALCON**

**AXYGEN**

**GOSELIN**

**PYREX**