

# Gyrozen 416 Benchtop Centrifuge Package

**VALID UNTIL:**  
31 March 2019

**Ref. Code:**  
GYRGZ-416

**The 416 Benchtop Centrifuge Package allows processing of sample tubes from 1.5 to 100ml using a 4-position swing-out rotor with 4x100ml buckets with sealing caps.**

The versatile Gyrozen 416 Centrifuge may be compact, but with 4x100ml loading capacity, it can accommodate a large variety of rotors, buckets, carriers and adaptors to meet your requirements. The unit is quiet in operation and safety is enhanced with imbalance cut out and lid lock protection.

- Max. Speed: 4,000rpm; Max Force: 2,774xg
- 10 program memory with selectable "Soft Start"
- Digital display of speed/force and time
- Many options and accessories including fixed angle/swing-out rotors + microplate carriers
- 2-year warranty

**A compact flexible low speed centrifuge for routine separations.**

**Speak to your local Bio-Strategy Representative for further details or call Customer Service on 0800 34 24 or email [sales.nz@bio-strategy.com](mailto:sales.nz@bio-strategy.com)**

Description	Code	List Price (excl. GST)	Sale Price (excl. GST)
416 Centrifuge + Rotor + 4 buckets + caps + adaptors	GYROGZ-416-Package	<b>\$5,405.00</b>	<b>\$4,575.00</b>

*Please note: Price includes centrifuge with swing-out rotor, 4x100ml buckets with sealing caps and 1 set of adaptors of your choice. Contact Bio-Strategy for adaptor options.*



Bio-Strategy Ltd.

T: 0800 34 24 66  
E: [sales.nz@bio-strategy.com](mailto:sales.nz@bio-strategy.com)  
[www.bio-strategy.com](http://www.bio-strategy.com)  
<https://shop.bio-strategy.com>

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.