

# Receive a 25% discount on a new SpectraMax ABS or ABS Plus Reader

VALID UNTIL:  
31 March 2019

Ref. Code:  
MDEV-ABS

## Affordable monochromator-based plate readers

Run various assays within the Vis/UV-Vis range such as ELISA, Bradford, BCA, microbial growth, protein quantitation and more with the new small but mighty ABS and ABS Plus Absorbance readers monochromator-based plate readers.

- ABS reader works with 96- and 384-well plates in the VIS range (340-850nm)
- ABS Plus works with 96-, 384 well plates as well as cuvettes in the UV/VIS range (190-1000nm)

### BENEFITS of the ABS and ABS PLUS Readers

- Flexibility to choose your wavelength, unlike filter-based readers
- Elimination of well-to-well pipetting differences with Patented PathCheck® Sensor
- Preconfigured protocols such as ELISA and Bradford in SoftMax® Pro Software
- Easy and fast assay setup and data analysis

Speak to your local Bio-Strategy Representative for further details or call Customer Service on 0800 34 24 or email [sales.nz@bio-strategy.com](mailto:sales.nz@bio-strategy.com)



**25% OFF,  
ORDER  
TODAY!**

Bio-Strategy Ltd.

T: 0800 34 24 66

E: [sales.nz@bio-strategy.com](mailto:sales.nz@bio-strategy.com)

[www.bio-strategy.com](http://www.bio-strategy.com)

<https://shop.bio-strategy.com>

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.