

# Pall Centrifugal Devices

## Save on Spin Devices

These devices are used for simple, reliable concentrating and desalting of samples.

Select the low protein binding Omega™ membrane for ultrafiltration applications. The Bio-Inert® membrane is ideal for microfiltration of biological solutions. The GHP membrane is ideal for sample preparation in analytical chemistry (HPLC, GC, IC and MS). It is suitable for aqueous, organic and aggressive organic samples and for rapid processing of samples.

- Rapid processing of samples
- Typical recoveries are greater than 90%
- Housing: Low binding PP
- Ultrasonically welded seals prevent bypass or seal failure

**Speak to your local Bio-Strategy Representative for further details or call Customer Service on 0800 34 24 or email [sales.nz@bio-strategy.com](mailto:sales.nz@bio-strategy.com)**

**10%  
OFF**

**VALID UNTIL:  
31 March 2019**

**Ref. Code:  
PALL-CENT**



Description	Pack Size	Code	Sale Price (excl. GST)
Nanosep® with Omega™ membrane, 3 kD, Grey	24 items	PALLOD003C33	<b>\$117.90</b>
Nanosep® with Omega™ membrane, 3 kD, Grey	100 items	PALLOD003C34	<b>\$393.30</b>
Nanosep® with Omega™ membrane, 3 kD, Grey	500 items	PALLOD003C35	<b>\$1,643.40</b>
Nanosep® with Omega™ membrane, 10 kD, Blue	24 items	PALLOD010C33	<b>\$117.90</b>
Nanosep® with Omega™ membrane, 10 kD, Blue	100 items	PALLOD010C34	<b>\$393.30</b>
Nanosep® with Omega™ membrane, 10 kD, Blue	500 items	PALLOD010C35	<b>\$1,647.00</b>
Nanosep® with Omega™ membrane, 30 kD, Red	24 items	PALLOD030C33	<b>\$116.10</b>
Nanosep® with Omega™ membrane, 30 kD, Red	100 items	PALLOD030C34	<b>\$392.32</b>
Nanosep® with Omega™ membrane, 30 kD, Red	500 items	PALLOD030C35	<b>\$1,643.40</b>
Nanosep® with Omega™ membrane, 100 kD, Clear	24 items	PALLOD100C33	<b>\$117.00</b>
Nanosep® with Omega™ membrane, 100 kD, Clear	100 items	PALLOD100C34	<b>\$393.30</b>
Nanosep® with Omega™ membrane, 100 kD, Clear	500 items	PALLOD100C35	<b>\$1,643.40</b>
Nanosep® with Omega™ membrane, 300 kD, Orange	24 items	PALLOD300C33	<b>\$117.00</b>
Nanosep® with Omega™ membrane, 300 kD, Orange	100 items	PALLOD300C34	<b>\$393.30</b>
Nanosep® with Omega™ membrane, 300 kD, Orange	500 items	PALLOD300C35	<b>\$1,643.40</b>
Nanosep® MF with GHP membrane, 0,45 µm, Clear	100 items	PALLODGHPC34	<b>\$414.90</b>
Nanosep® MF with GHP membrane, 0,45 µm, Clear	500 items	PALLODGHPC35	<b>\$1,410.30</b>
Nanosep® MF with Bio-Inert® membrane, 0,2 µm, Turquoise	24 items	PALLODM02C33	<b>\$86.40</b>
Nanosep® MF with Bio-Inert® membrane, 0,2 µm, Turquoise	100 items	PALLODM02C34	<b>\$300.60</b>
Nanosep® MF with Bio-Inert® membrane, 0,2 µm, Turquoise	500 items	PALLODM02C35	<b>\$1,286.10</b>
Nanosep® MF with Bio-Inert® membrane, 0,45 µm, Purple	24 items	PALLODM45C33	<b>\$86.40</b>
Nanosep® MF with Bio-Inert® membrane, 0,45 µm, Purple	100 items	PALLODM45C34	<b>\$300.60</b>
Nanosep® MF with Bio-Inert® membrane, 0,45 µm, Purple	500 items	PALLODM45C35	<b>\$1,286.10</b>

Bio-Strategy Ltd.

T: 0800 34 24 66

E: [sales.nz@bio-strategy.com](mailto:sales.nz@bio-strategy.com)

[www.bio-strategy.com](http://www.bio-strategy.com)

<https://shop.bio-strategy.com>

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.