

## Next Generation CO<sub>2</sub> Incubators Providing Optimum Cell Culture Conditions

### Grow Results – Not Bacteria

**PHCbi CO<sub>2</sub> incubators combine both background and active features to ensure you're not compromising your research.**

inCu-saFe is PHCbi's background contamination control. A copper-enriched stainless-steel alloy gives you the best of both worlds; the germicidal properties of copper without the gradual oxidation-discoloration, to keep your interior looking like-new for a longer time!

Choose from either 165L, 170L or 230L models. Base unit (MC-18AC) features proven reliable TC sensor whilst the upgradeable models have the ceramic based single beam dual wavelength IR sensor with continuous auto-zero calibration!

- inCu-saFe – copper enriched stainless-steel alloy interior and shelves
- DHA – direct heat & air jacket with 3-zone heating for precise temperature control
- CO<sub>2</sub> range – 0 to 20%
- Temperature range – ambient + 5 to 50°C
- Field reversible door

**Speak to your local Bio-Strategy Representative for further details or call Customer Service on 0800 34 24 or email [sales.nz@bio-strategy.com](mailto:sales.nz@bio-strategy.com)**

Description	Code	List Price (excl. GST)	Sale Price (excl. GST)
CO <sub>2</sub> Incubator (170L) base unit with TC sensor	SANYMCO-18AC	<b>\$5,825.00</b>	<b>\$4,995.00</b>
CO <sub>2</sub> Incubator (165L) base unit with IR sensor (upgradeable)	SANYMCO-170AICL-PE	<b>\$7,118.00</b>	<b>\$6,265.00</b>
CO <sub>2</sub> Incubator (165L) base unit with IR sensor, SafeCell UV decontamination + H <sub>2</sub> O <sub>2</sub> decon board + door lock (upgradeable)	SANYMCO-170AICUVHL-PE	<b>\$9,265.00</b>	<b>\$8,155.00</b>
CO <sub>2</sub> Incubator (230L) base unit with IR sensor, SafeCell UV decontamination (upgradeable)	SANYMCO-230AICUVL-PE	<b>\$8,925.00</b>	<b>\$7,855.00</b>



**VALID UNTIL:  
31 March 2019**

**Ref. Code:  
SANY-CO2**

**In Stock &  
12% Off!**

Bio-Strategy Ltd.

T: 0800 34 24 66  
E: [sales.nz@bio-strategy.com](mailto:sales.nz@bio-strategy.com)  
[www.bio-strategy.com](http://www.bio-strategy.com)  
<https://shop.bio-strategy.com>

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.