

Save 30% on QIAseq FastSelect Kits

One-step rRNA removal in just 10 seconds and without bead cleanup? It's possible with QIAseq FastSelect!

1 step vs. 37 steps – take a shortcut to RNA-seq success with QIAseq FastSelect

New QIAseq FastSelect RNA Removal Kits are the latest revolution rocking the world of RNA-seq, providing the fastest ever removal of unwanted RNA from your libraries! Why spend hours removing ribosomal RNA (rRNA) and globin mRNA from your RNA-seq library when you can achieve it in **just 20 minutes using a single pipetting step during library prep – and with only 10 seconds of hands-on time!**

Sequence what matters – utilise your sequencing resources and capacity more efficiently

Even if you're working with fragmented RNA samples, our radically simplified workflow delivers exceptional quality and performance compared to traditional techniques such as magnetic bead hybridization and capture or enzymatic methods. This means you can speed up RNA-seq experiments and significantly reduce read allocation to unwanted RNAs – while attaining reproducible RNA sequencing results.

Compatible with QIAGEN, Illumina, NEB and KAPA RNA stranded library kits.

QIAseq FastSelect Custom RNA Removal Kits designed for your RNA species of choice is available on request.

VALID UNTIL:

31 March 2019

Ref. Code:

FastSelect19

Rethink
RNA-seq.

**30%
OFF!**

Description	Pack Size	Code	List Price (excl. GST)	Sale Price (excl. GST)
QIAseq FastSelect RNA Removal kit-Human,Cyto & Mito rRNA,24	24	QIAG333180-THS-001Z-24	\$2,834.76	\$1,984.30
QIAseq FastSelect RNA Removal kit-Human,Cyto & Mito rRNA,96	96	QIAG333180-THS-001Z-96	\$10,196.47	\$7,137.41
QIAseq FastSelect RNA Removal kit-Mouse,Cyto & Mito rRNA,24	24	QIAG333180-TMM-001Z-24	\$2,834.76	\$1,984.30
QIAseq FastSelect RNA Removal kit-Mouse,Cyto & Mito rRNA,96	96	QIAG333180-TMM-001Z-96	\$10,196.47	\$7,137.41
QIAseq FastSelect RNA Removal kit-Rat,Cyto & Mito rRNA,24	24	QIAG333180-TRN-001Z-24	\$2,834.76	\$1,984.30
QIAseq FastSelect RNA Removal kit-Rat,Cyto & Mito rRNA,96	96	QIAG333180-TRN-001Z-96	\$10,196.47	\$7,137.41
QIAseq FastSelect RNA Removal kit-Human,Globin RNA,24	24	QIAG333180-THS-002Z-24	\$2,834.76	\$1,984.30
QIAseq FastSelect RNA Removal kit-Human,Globin RNA,96	96	QIAG333180-THS-002Z-96	\$10,196.47	\$7,137.41
QIAseq FastSelect RNA Removal kit-Mouse,Globin RNA,24	24	QIAG333180-TMM-002Z-24	\$2,834.76	\$1,984.30
QIAseq FastSelect RNA Removal kit-Mouse,Globin RNA,96	96	QIAG333180-TMM-002Z-96	\$10,196.47	\$7,137.41
QIAseq FastSelect RNA Removal kit-Rat,Globin RNA,24	24	QIAG333180-TRN-002Z-24	\$2,834.76	\$1,984.30
QIAseq FastSelect RNA Removal kit-Rat,Globin RNA,96	96	QIAG333180-TRN-002Z-96	\$10,196.47	\$7,137.41
QIAseq FastSelect Multi-RNA Removal,Human,rRNA&GlobinmRNA,24	24	QIAG333280-THS-201Z-24	\$3,883.34	\$2,718.29
QIAseq FastSelect Multi-RNA Removal,Human,rRNA&GlobinmRNA,96	96	QIAG333280-THS-201Z-96	\$14,642.41	\$10,249.53
QIAseq FastSelect Multi-RNA Removal,Human,rRNA&GlobinmRNA,24	24	QIAG333280-TMM-201Z-24	\$3,883.34	\$2,718.29
QIAseq FastSelect Multi-RNA Removal,Human,rRNA&GlobinmRNA,96	96	QIAG333280-TMM-201Z-96	\$14,642.41	\$10,249.53
QIAseq FastSelect Multi-RNA Removal,Human,rRNA&GlobinmRNA,24	24	QIAG333280-TRN-201Z-24	\$3,883.34	\$2,718.29
QIAseq FastSelect Multi-RNA Removal,Human,rRNA&GlobinmRNA,96	96	QIAG333280-TRN-201Z-96	\$14,642.41	\$10,249.53
QIAseq Stranded mRNA Select Kit (24)	24	QIAG180773	\$2,244.67	\$1,571.24
QIAseq Stranded mRNA Select Kit (96)	96	QIAG180775	\$7,138.97	\$4,997.20
QIAseq Stranded Total Lib Kit (24)	24	QIAG180743	\$2,120.29	\$1,484.18
QIAseq Stranded Total Lib Kit (96)	96	QIAG180745	\$6,888.76	\$4,822.06

Bio-Strategy Ltd.

T: 0800 34 24 66

E: sales.nz@bio-strategy.com

www.bio-strategy.com

<https://shop.bio-strategy.com>

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.