

## The Most Economic VWR Pipettor EXCLUSIVE TO BIO-STRATEGY

### Time to replace your old pipettors!

#### VWR Standard Line Economic Pipettors

These lightweight Standard Line pipettes provide accuracy and precision. Ergonomic features include contoured handle shape and a soft spring system. 3 years warranty, autoclavable.

- Easily adjustable ejector height system for many tips
- UV resistant and easy to calibrate
- Colour coded shaft for easy volume identification

Description	Accuracy	Imprecision	Availability	Code	Sale Price (excl. GST)
<b>Single channel pipette, mechanical, Standard Line</b>					
0.2-2µL	±12.0-±1.5%	≤6.0-≤0.7%	stocked	VWRI613-2891	<b>\$207.00</b>
0.5-10µL	±4.0-±0.5%	≤2.8-≤0.4%	stocked	VWRI613-2892	<b>\$207.00</b>
2-20µL	±3.0-±0.8%	≤1.5-≤0.3%	enquire	VWRI613-2893	<b>\$207.00</b>
10-100µL	±1.6-±0.8%	≤0.8-≤0.2%	stocked	VWRI613-2894	<b>\$207.00</b>
20-200µL	±1.2-±0.6%	≤0.6-≤0.2%	stocked	VWRI613-2895	<b>\$207.00</b>
100-1000µL	±1.6-±0.6%	≤0.4-≤0.15%	stocked	VWRI613-2896	<b>\$207.00</b>
500-5000µL	±1.2-±0.5%	≤0.5-≤0.15%	stocked	VWRI613-2897	<b>\$214.00</b>
<b>Packs</b>					
WR Standard 3 Pack (2-20, 20-200 and 100-1000µl)			enquire	VWRI613-5469	<b>\$603.00</b>

#### Bundle up with VWR Next Generation Pipette Tips for further savings!



**VALID UNTIL:**  
31 March 2019  
**Ref. Code:**  
VWRI-PIPE

**Only  
\$207 per  
pipettor!**



Bio-Strategy Ltd.

T: 0800 34 24 66  
E: [sales.nz@bio-strategy.com](mailto:sales.nz@bio-strategy.com)  
[www.bio-strategy.com](http://www.bio-strategy.com)  
<https://shop.bio-strategy.com>

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.