

Falcon FACS Tubes and Multi-Flasks

VALID UNTIL:
30 June 2019

Ref. Code:
BDAA-CELL

Falcon® Round Bottom Tubes

Designed for flow cytometry applications, the Falcon® 12x75mm tube with cell strainer cap offers a convenient way to prepare laboratory samples. A 35µm nylon mesh is incorporated into the dual-position snap cap, which can be used to collect the dissociated sample for downstream, processing in flow cytometry instruments. RNase/DNase-free.

- Widely used and referenced in laboratory protocols
- Dual position snap-cap prevents sample loss with ease of handling
- Polystyrene tubes are best suited for procedures requiring high optical clarity

Falcon® Cell Culture Multi-Flasks

Increase productivity by enabling to grow more cells faster and easier. 3- or 5-layer formats provide 525cm² and 875cm² cell growth surface area, respectively.

- Even distribution of media across all layers for homogeneous cell growth
- Ability to mix cells and reagents in the Falcon Multi-Flask saves time and reduces risk of contamination
- Flexible design allows to pour or aspirate/recover cells using a pipette
- Consistent surface treatment for predictable scale-up
- Lot number printed on each flask for traceability

Speak to your local Bio-Strategy Representative for further details or call Customer Service on 0800 34 24 66 or email sales.nz@bio-strategy.com



**Falcon®
Cell Culture
Specials**



Description	Pack Size	Code	Sale Price (excl. GST)
Falcon® 5ml Round Bottom Polystyrene Test Tube, with Cell Strainer Snap Cap	500	BDAA352235	\$542.00
Falcon® 5ml Round Bottom High Clarity PolyPropylene Test Tube, with Snap Cap, Sterile	500	BDAA352063	\$200.00
Falcon® 5ml Round Bottom PolyPropylene Test Tube, without Cap, Sterile	1,000	BDAA352053	\$184.00
Falcon® 5ml Round Bottom PolyPropylene Test Tube, without Cap	1,000	BDAA352002	\$105.00
Falcon® 5ml Round Bottom Polystyrene Test Tube, with Snap Cap, Sterile, Individually Wrapped	500	BDAA352003	\$181.00
Falcon® 5ml Round Bottom Polystyrene Test Tube, with Snap Cap, Sterile	500	BDAA352058	\$174.00
Falcon® 5ml Round Bottom Polystyrene Test Tube, with Snap Cap, Sterile,	1,000	BDAA352054	\$153.00
Falcon® 5ml Round Bottom Polystyrene Test Tube, without Cap, Sterile	1,000	BDAA352052	\$231.00
Falcon® 5ml Round Bottom Polystyrene Test Tube, without Cap, Non-Sterile	1,000	BDAA352008	\$104.00
Falcon® Snap Caps for 12x75mm Test Tubes, Sterile	2,000	BDAA352032	\$159.00
Falcon® 525cm ² Cell Culture Multi-Flask, 3-Layer with Vented Cap	12	BDAA353143	\$143.00
Falcon® 875cm ² Cell Culture Multi-Flask, 5-Layer with Vented Cap	8	BDAA353144	\$172.00

Bio-Strategy Ltd.

T: 0800 34 24 66

E: sales.nz@bio-strategy.com

www.bio-strategy.com

<https://shop.bio-strategy.com>

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.