


## Bel-Art No-Wire Test Tube Racks

### Submersible! Stackable! Colour-Coded! Autoclavable!

A great alternative to wire racks, with price and space saving advantages over other plastic racks. Dimensions for all racks are 24.6x10.4x6.4cm high (9.7x4.1x2.5") and they assemble in seconds.

- Made of a special blend of autoclavable polypropylene; will sink in a water bath and maintain their stability even in agitated water
- Also perfect for use outside the water bath; meets the demands of daily usage in the laboratory
- Available in numerous hole sizes, with a choice of three colours for safety and identification
- Features end plates for labelling; numbered and lettered rows for tube identification
- Interlocking feet allow stacking of racks for storage
- Steam autoclavable at 121°C (250°F)

**Speak to your local Bio-Strategy Representative for further details or call Customer Service on 0800 34 24 66 or email [sales.nz@bio-strategy.com](mailto:sales.nz@bio-strategy.com)**



**Special  
Sale  
Prices!**

Description	Pack Size	Code	Sale Price (excl. GST)
For 10-13mm Tubes, 90 Place, White	1 Item	BELAF187450000	<b>\$24.00</b>
For 10-13mm Tubes, 90 Place, Red	1 Item	BELAF187460000	<b>\$24.00</b>
For 10-13mm Tubes, 90 Place, Blue	1 Item	BELAF187470000	<b>\$24.00</b>
For 10-13mm Tubes, 60 Place, White	1 Item	BELAF187450001	<b>\$24.00</b>
For 13-16mm Tubes, 60 Place, Red	1 Item	BELAF187460001	<b>\$24.00</b>
For 13-16mm Tubes, 60 Place, Blue	1 Item	BELAF187470001	<b>\$24.00</b>
For 13-16mm Tubes, 40 Place, White	1 Item	BELAF187450002	<b>\$24.00</b>
For 16-20mm Tubes, 40 Place, Red	1 Item	BELAF187460002	<b>\$24.00</b>
For 16-20mm Tubes, 40 Place, Blue	1 Item	BELAF187470002	<b>\$24.00</b>
For 16-20mm Tubes, 24 Place, White	1 Item	BELAF187450003	<b>\$24.00</b>
For 20-25mm Tubes, 24 Place, Red	1 Item	BELAF187460003	<b>\$24.00</b>
For 20-25mm Tubes, 24 Place, Blue	1 Item	BELAF187470003	<b>\$24.00</b>
For 20-25mm Tubes, 21 Place, White	1 Item	BELAF187450004	<b>\$24.00</b>
For 25-30mm Tubes, 21 Place, Red	1 Item	BELAF187460004	<b>\$24.00</b>
For 25-30mm Tubes, 21 Place, Blue	1 Item	BELAF187470004	<b>\$24.00</b>

Bio-Strategy Ltd.

T: 0800 34 24 66

E: [sales.nz@bio-strategy.com](mailto:sales.nz@bio-strategy.com)

[www.bio-strategy.com](http://www.bio-strategy.com)

<https://shop.bio-strategy.com>

**VALID UNTIL:  
30 June 2019**

**Ref. Code:  
BELA-RACKS**

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.