

Quickfit® Precision Joints, Outstanding Performance - 25% OFF!

Top Quality Interchangeable Jointed Glassware on Sale Now!

All Quickfit® glass items are exclusively manufactured from Pyrex® borosilicate glass. The Quickfit® range comprises of around 1,000 individual products that can be combined to construct an extensive range of assemblies.

- Low coefficient of expansion for use over a wide temperature range
- Highly resistant to chemical attack from most acids, alkalis and organic solvents
- 2-stage grinding process for an exceptional seal and reduced risk of joint seizure
- Excellent surface finish to prevent leaking and ensure gas and vacuum tightness
- Rounded rim profile for extra joint strength and reduced chipping
- 25% off list on all Quickfit® Interchangeable Jointed Glassware!

Speak to your local Bio-Strategy Representative for further details or call Customer Service on 0800 34 24 66 or email sales.nz@bio-strategy.com



**25% off
Quickfit®
Glassware**

Bio-Strategy Ltd.

T: 0800 34 24 66

E: sales.nz@bio-strategy.com

www.bio-strategy.com

<https://shop.bio-strategy.com>

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.

**VALID UNTIL:
30 June 2019**

**Ref. Code:
QUICKFIT**