

ESCO Laminar Flow Benches

Choose the leader in cleanroom equipment solutions!

The new generation Airstream Laminar Flow Clean Benches offer a wide range of options; with glass side walls, sizes ranging from 3 to 8 feet and an internal height from 2 feet, 2.5 feet or 3 feet. The full range of the new Esco Airstream Laminar Flow Clean Benches are powered by the latest generation DC ECM motor. They are more efficient than the legacy ECM and VFD motors. They now also offer ISO Class 3 workzone with 10x efficiency of its HEPA filters. With more choices, the Esco Laminar Flow Clean Benches still offer the same level of product protection for your samples and processes where operator protection is not required.

Key Benefits:

- World leading Low Noise Laminar Flow Cabinet
- Energy-efficient DC ECM motor with power consumption of only 155 watts
- With standby mode for more energy saving
- Stable and self-compensating flow
- With ISOCIDETM antibacterial powder coating
- ULPA filter (up to ISO Class 3)

Speak to your local Bio-Strategy Representative for further details or call Customer Service on 0800 34 24 66 or email sales.nz@bio-strategy.com



**SAVE
10%**

**VALID UNTIL:
30 June 2019**

**Ref. Code:
ESCO-FLOW**

Bio-Strategy Ltd.

T: 0800 34 24 66

E: sales.nz@bio-strategy.com

www.bio-strategy.com

<https://shop.bio-strategy.com>

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.