

Hamilton Company Neuros Syringes

VALID UNTIL:
30 June 2019

Ref. Code:
HAMC-NEUR

Maximum Control, Minimum Damage

Make sure your syringe is the finest precision fluid measuring instrument in the world.

- Neuros syringe technology provides unprecedented functionality for controlled animal injections. An industry first, the Neuros accurately dispenses 50nl to 100µl of liquid through an ultrafine needle protected by a custom-designed sleeve. Developed specifically for the neurosciences, the Neuros minimises injection site damage by increasing needle rigidity and delivering a microvolume of liquid to an exact location.
- Needle rigidity improves insertion path accuracy
- Minimal tissue damage reduces injection variability
- Reduced sample loss saves money and materials
- Fine gauge needle creates smaller injection sites
- Compatibility with most infusion pumps and stereotaxic holders means an easy integration into existing processes

**20%
OFF**



Speak to your local Bio-Strategy Representative for further details or call Customer Service on 0800 34 24 66 or email sales.nz@bio-strategy.com

Description	Pack Size	Code	List Price (excl. GST)	Sale Price (excl. GST)
0.5µl, Neuros Model 7000.5 KH Syringe, Point Style 3	Each	HAMC65457-01	\$600.46	\$480.37
0.5µl, Neuros Model 7000.5 KH Syringe, Point Style 4	Each	HAMC65457-02	\$600.46	\$480.37
1.0µl, Neuros Model 7001 KH Syringe, Point Style 3	Each	HAMC65458-01	\$600.46	\$480.37
1.0µl, Neuros Model 7001 KH Syringe, Point Style 4	Each	HAMC65458-02	\$600.46	\$480.37
2.0µl, Neuros Model 7002 KH Syringe, Point Style 3	Each	HAMC65459-01	\$600.46	\$480.37
2.0µl, Neuros Model 7002 KH Syringe, Point Style 4	Each	HAMC65459-02	\$600.46	\$480.37
5µl, Neuros Model 75 RN Syringe, Point Style 3	Each	HAMC65460-02	\$516.40	\$413.12
5µl, Neuros Model 75 RN Syringe, Point Style 4	Each	HAMC65460-03	\$516.40	\$413.12
10µl, Neuros Model 1701 RN Syringe, Point Style 3	Each	HAMC65460-05	\$480.37	\$384.29
10µl, Neuros Model 1701 RN Syringe, Point Style 4	Each	HAMC65460-06	\$480.37	\$384.29
25µl, Neuros Model 1702 RN Syringe, Point Style 3	Each	HAMC65460-10	\$480.37	\$384.29
25µl, Neuros Model 1702 RN Syringe, Point Style 4	Each	HAMC65460-11	\$480.37	\$384.29
50µl, Neuros Model 1705 RN Syringe, Point Style 3	Each	HAMC65460-15	\$480.37	\$384.29
50µl, Neuros Model 1705 RN Syringe, Point Style 4	Each	HAMC65460-16	\$480.37	\$384.29
100µl, Neuros Model 1710 RN Syringe, Point Style 3	Each	HAMC65460-20	\$480.37	\$384.29
100µl, Neuros Model 1710 RN Syringe, Point Style 4	Each	HAMC65460-21	\$480.37	\$384.29

Bio-Strategy Ltd.

T: 0800 34 24 66
E: sales.nz@bio-strategy.com
www.bio-strategy.com
<https://shop.bio-strategy.com>

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.