

Hamilton Plate Sealer - Standalone or Integrate with your Hamilton Robot

Purchase a Hamilton Plate Sealer and receive EITHER a free integration kit OR \$3,240.00 off the list price!

Compact. Flexible. Economical.

The Hamilton Plate Sealer is designed to seal standard ANSI/SLAS footprint plates with optically clear or pierceable foil seals. The Plate Sealer's flexible design allows it to seal plates 8–45mm in height without using compressed air or thermal control. Standalone or integration – the choice is yours!

As a standalone device, the Plate Sealer offers efficient and reliable placement of seals with the convenience of an intuitive touch screen. The Plate Sealer automates the sealing process, while the operator loads plates and seals, controls the device using touch screen prompts, and unloads plates once sealed.

Integration offers the convenience of additional walk-away time as plates and seals are automatically supplied to the seamlessly integrated on-deck sealer.

Speak to your local Bio-Strategy Representative for further details or call Customer Service on 0800 34 24 66 or email sales.nz@bio-strategy.com

Description	Code	List Price (excl. GST)	Sale Price (excl. GST)
Hamilton Plate Sealer	HAMR93005-01	\$25,920.00	\$22,680.00*
Standard Integration Kit, STAR/Vantage	HAMR97553-01	\$3,240.00	No Charge**
Aluminium Seal, 2x50/Pack	HAMR67760-01	\$235.00	\$216.00
Optically Clear Seal, Non-Pierceable, 2x50/Pack	HAMR67765-01	\$496.00	\$470.00

* Offer 1: \$3,240.00 off list price for a standalone purchase

** Offer 2: Free Standard Integration Kit valued at \$3,240.00, available for owners of a Hamilton Robot. Fitting service if required is at an extra cost. Please inquire about details specific to your Hamilton Robot model.



Bio-Strategy Ltd.

T: 0800 34 24 66

E: sales.nz@bio-strategy.com

www.bio-strategy.com

<https://shop.bio-strategy.com>

VALID UNTIL:
30 June 2019

Ref. Code:
HAMRSEAL

**Free
Standard
Integration
Kit or \$3,240
OFF!**

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.