

# Heidolph Hei-Tec Magnetic Stirrers on Sale

**Safe heating and mixing with uncompromising quality and simple handling.**

The robust magnetic stirrers of the Hei-PLATE series were developed to optimally mix low-viscosity liquids from smooth to intense. The premium devices offer the highest operating safety and comfort standards. 800W (for 230V models) heating capacity and the special Kera-Disk® hotplate with aluminium core allow for minimum heat-up phases that result in lasting cost reductions.

- This unit features a digital display and allows for easy setting and monitoring of values in the illuminated readout
- For your protection an independent safety circuit will switch off heating if the hotplate temperature exceeds 25°C over the set temperature
- Digital speed setting from 100 to 1,400rpm at an accuracy of  $\pm 2\%$  and temperature setting of up to 300°C

A separate on/off switch for heating prevents unintentional heat-up. If heating is switched on the on/off button is illuminated for visual control. In addition, a residual heat indicator prevents from burning hazards when the heat function is switched off. An illuminated button indicates clearly if the stirring mode is activated.

Upgrade this magnetic stirrer with the optional Pt1000 temperature sensor for precise temperature control, overshoot protection and reproducible results.

**Speak to your local Bio-Strategy Representative for further details or call Customer Service on 0800 34 24 66 or email [sales.nz@bio-strategy.com](mailto:sales.nz@bio-strategy.com)**



**20% OFF  
STOCKED  
UNITS**

Bio-Strategy Ltd.

T: 0800 34 24 66

E: [sales.nz@bio-strategy.com](mailto:sales.nz@bio-strategy.com)

[www.bio-strategy.com](http://www.bio-strategy.com)

<https://shop.bio-strategy.com>

**VALID UNTIL:  
30 June 2019**

**Ref. Code:  
HEI-TEC**

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.