

Heidolph Hei-Torque Core Overhead Stirrers on Sale

Lightweight for big tasks.

The compact design allows for integration in closed systems, such as fume hoods, reactors, or production systems. Suitable for low- to medium-viscosity media up to 25L. The large diameter of the chuck (10mm) allows you to use even larger impellers and VISCO JET® tools. This facilitates a wide variety of applications, such as homogenisation, dispersing, the dissolving of agglomerates and more.

- Speed range up to 2,000rpm
- Viscosity up to 10,000mPas
- Control knob for rotation speed, pushing starts or stops the function
- Timer
- “Max” button for short-term operation at maximum speed

Speak to your local Bio-Strategy Representative for further details or call Customer Service on 0800 34 24 66 or email sales.nz@bio-strategy.com

**VALID UNTIL:
30 June 2019**

**Ref. Code:
HEI-TORQUE**

**20% OFF
STOCKED
UNITS**



Bio-Strategy Ltd.

T: 0800 34 24 66

E: sales.nz@bio-strategy.com

www.bio-strategy.com

<https://shop.bio-strategy.com>

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.