

# Kimble Plain Disposable Borosilicate Culture Tubes on Sale

## Premium Quality Borosilicate Culture Tubes on Sale – Order Today!

Kimble Plain Disposable Premium Quality Borosilicate Glass Culture Tubes with sturdy, uniform bottoms and consistent lengths.

- 51 Expansion borosilicate glass for excellent resistance against chemical leaching
- Conforms to USP Type I and ASTM E438, Type I, Class B requirements
- Designed from ASTM Specification E890
- Packed in shrink wrapped trays for ease of use and protection from breaking and contamination
- No marking spots
- Polypropylene colour-coded, autoclavable KIM-KAP™ closures available


**Speak to your local Bio-Strategy Representative for further details or call Customer Service on 0800 34 24 66 or email [sales.nz@bio-strategy.com](mailto:sales.nz@bio-strategy.com)**

Description	Pack Size	Availability	Code	List Price (excl. GST)	Sale Price (excl. GST)
Tube Plain Disposable Culture Borosilicate 6x50mm (known as "Durham Tubes")	100 Items	Enquire	KIML73500-650	<b>\$13.98</b>	<b>\$9.79</b>
Tube Plain Disposable Culture Borosilicate 10x75mm, Suitable for Cell Washing Procedures	1,000 Items	In Stock	KIML73500-1075	<b>\$113.71</b>	<b>\$79.60</b>
Tube Plain Disposable Culture Borosilicate 12x75mm, Suitable for Cell Washing Procedures	1,000 Items	In Stock	KIML73500-1275	<b>\$117.73</b>	<b>\$82.41</b>
Tube Plain Disposable Culture Borosilicate 13x100mm	1,000 Items	Enquire	KIML73500-13100	<b>\$147.71</b>	<b>\$103.40</b>
Tube Plain Disposable Culture Borosilicate 16x100mm	1,000 Items	In Stock	KIML73500-16100	<b>\$197.91</b>	<b>\$138.53</b>
Tube Plain Disposable Culture Borosilicate 16x125mm (formerly known as 60K18B Tubes)	1,000 Items	In Stock	KIML73500-16125	<b>\$230.18</b>	<b>\$161.12</b>
Tube Plain Disposable Culture Borosilicate 16x150mm	1,000 Items	In Stock	KIML73500-16150	<b>\$248.14</b>	<b>\$173.70</b>
Tube Plain Disposable Culture Borosilicate 18x150mm (formerly known as 60F165B Tubes)	500 Items	Enquire	KIML73500-18150	<b>\$185.05</b>	<b>\$129.53</b>
Tube Plain Disposable Culture Borosilicate 20x150mm	500 Items	Enquire	KIML73500-20150	<b>\$268.04</b>	<b>\$187.63</b>
Tube Plain Disposable Culture Borosilicate 25x150mm	500 Items	Enquire	KIML73500-25150	<b>\$449.64</b>	<b>\$314.75</b>

**VALID UNTIL:  
30 June 2019**

**Ref. Code:  
KIML-TUBES**

**Buy Now  
and Save!**



Bio-Strategy Ltd.

T: 0800 34 24 66

E: [sales.nz@bio-strategy.com](mailto:sales.nz@bio-strategy.com)

[www.bio-strategy.com](http://www.bio-strategy.com)

<https://shop.bio-strategy.com>

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.