

# Starna Scientific

## The Spectroscopy Specialists

### Cell/Cuvettes for all Spectrophotometer, Fluorimeter, Polarimeter and Laser applications.

Starna Scientific offer a wide range of Cell Sizes from Standard Macro Cells down to Ultra-Micro Cells and Flow Cells. Choose from a range of Glass and Quartz Materials to cover the spectrum from 190nm (Far UV) through to 3800nm (Near Infra-Red). Starna also offer a range of NIST Traceable Glass and Liquid reference materials to meet your Ph. Eur., USP, BP requirements.

- Starna Cells fused method of construction, ensures maximum physical strength as well as resistance to solvents
- Most cells can be used safely with pressure differentials of up to  $3 \times 10^5$  Pa and some up to  $10 \times 10^5$  Pa
- Unconditional guarantee against faulty workmanship
- Matched cells available
- 10% discount off list price for all Starna cells – ask for a quote today

**Speak to your local Bio-Strategy Representative for further details or call Customer Service on 0800 34 24 66 or email [sales.nz@bio-strategy.com](mailto:sales.nz@bio-strategy.com)**

Description	Pack Size	Availability	Code	List Price (excl. GST)	Sale Price (excl. GST)
Cell Type 1 Quartz Far UV 10mm (2 Match Cells/Box)	2 Items	In Stock	OPTI1/Q/10/M2	\$191.46	\$162.74



**10%  
OFF ALL  
Starna  
Cells**

Bio-Strategy Ltd.

T: 0800 34 24 66

E: [sales.nz@bio-strategy.com](mailto:sales.nz@bio-strategy.com)

[www.bio-strategy.com](http://www.bio-strategy.com)

<https://shop.bio-strategy.com>

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.

**VALID UNTIL:  
30 June 2019**

**Ref. Code:  
OPTI-CELLS**