

# Pall Multi-Well Plate Manifold and Vacuum Pump

**VALID UNTIL:**  
30 June 2019

**Ref. Code:**  
PALL-PLATE

**Designed to Perfectly Fit AcroPrep™, AcroPrep Advance, and AcroWell™ Filter Plates.**

## Optimise and simplify your plate vacuum applications!

Comes complete with the necessary O-ring and gasket. The control block includes the vacuum pressure gauge, vacuum metering valve, vacuum release valve, and 1/4 in. hose barb for vacuum line attachment.

Vacuum manifold unit includes a Delrins spacer block designed to accommodate standard 350µl receiver plates. The spacer block has been optimised to reduce the space between the receiver plate and the filter plate during vacuum filtration.

- Optional spacer blocks are available for use with 1ml receiver plates
- Adapter collar holds filter plates tightly to receiver plates for centrifugation
- The multi-well plate vacuum manifold is an anodised aluminium manifold that has been designed and optimised for the vacuum filtration of AcroPrep, AcroPrep Advance, and AcroWell multi-well filter plates

**Speak to your local Bio-Strategy Representative for further details or call Customer Service on 0800 34 24 66 or email [sales.nz@bio-strategy.com](mailto:sales.nz@bio-strategy.com)**

Description	Pack Size	Code	Sale Price (excl. GST)
Multi-Well Plate Manifold	1 Item	PALL5017	<b>\$1,644.00</b>
Vacuum/Pressure Pump	1 Item	PALL13158	<b>\$1,609.00</b>

**ON SALE!**



Bio-Strategy Ltd.

T: 0800 34 24 66

E: [sales.nz@bio-strategy.com](mailto:sales.nz@bio-strategy.com)

[www.bio-strategy.com](http://www.bio-strategy.com)

<https://shop.bio-strategy.com>

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.