

PHCbi MDF-U74V -86°C Upright Freezer, 728 Litres on Sale

Proven reliable -86°C storage with cascade refrigeration and VIP Insulation.

PHCbi VIP Freezers ensures maximum storage capacity with a compact footprint, combined with vacuum insulation panel technology.

The PHCbi MDF-U74V continues the legacy of Ultra-Low Temperature -86°C freezers from the former Sanyo/Panasonic stable with vacuum panel insulation technology (VIP) and Cool Safe Compressors to provide a reliable unit with increased compressor longevity and lower heat output. The units are quiet in operation, run off a standard 10-amp supply and the easy access airflow filter reduces maintenance times.

- Cascade refrigeration with **Cool Safe compressors:** proven reliable cooling technology with compressors running at lower temperatures
- **VIP Plus insulation:** provides up to 30% more storage capacity compared to conventionally insulated freezers with the same footprint
- User friendly eye-level LED display with audible and visual alarms
- Two insulated inner doors with eagle-latches reduces heat loss
- Quiet in operation: 51dB(A) and low power consumption: 8.5kWh/24h
- Keyed door lock for security

Speak to your local Bio-Strategy Representative for further details or call Customer Service on 0800 34 24 66 or email sales.nz@bio-strategy.com

Description	Code	Usual Sell (excl. GST)	Special Price (excl. GST)
-86°C Upright Freezer, 728 litres	SANYMDF-U74V	\$20,547.00	\$15,975.00

**VALID UNTIL:
30 June 2019**

**Ref. Code:
MDF-U74V**

**Special Price:
\$15,975.00**



Bio-Strategy Ltd.

T: 0800 34 24 66
E: sales.nz@bio-strategy.com
www.bio-strategy.com
<https://shop.bio-strategy.com>

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.