

Scilabware Azlon® Magnetic Stirrer and Spin Bars

VALID UNTIL:
30 September 2019

Ref. Code:
AZLO-BARS

Magnetic PTFE Stirrer Bars in a variety of shapes and sizes.

For the stirring of liquids on magnetic stirrers and hot plate stirrers. Strong alloy magnets sheathed in pure white inert PTFE. Supplied in polarity aligned packs to ensure maximum coupling strength. Resistant to virtually all chemical attack and temperatures up to 260°C. The retrievers have a magnet permanently embedded in the PTFE rod and are ideal for retrieving bars from particularly sensitive or corrosive liquids.

Choose the bar that's right for you or an assorted box of the most popular sizes and shapes

- Cylindrical plain with rounded ends, suitable for a wide range of flat-bottomed containers
- Plain similar to cylindrical bars but give greater turbulence at low speed, ideal for use with large vessels
- Oval for efficient stirring in round bottom flasks
- Central pivot ring ensures efficient spinning even on curved or uneven bases
- Octagonal with central pivot ring, the shape gives increased turbulence at low speeds
- Cross head design single or double sided, for efficient stirring in tube-like containers and effective as general stirrers
- Micro for use in small vessels, e.g. test tubes, vials etc.

Speak to your local Bio-Strategy Representative for further details or call Customer Service on 0800 34 24 66 or email sales.nz@bio-strategy.com

**On
Special,
Order
Today!**



Description	Pack Size	Availability	Code	Sale Price (excl. GST)
Stir Bar Pivot Ring 8x3mm PTFE	2 Items	Enquire	AZLOSWN500	\$4.01
Stir Bar Pivot Ring 15x4.5mm PTFE	2 Items	Enquire	AZLOSWN501	\$4.01
Stir Bar Pivot Ring 12x4.5mm PTFE	2 Items	Enquire	AZLOSWN502	\$3.76
Stir Bar Pivot Ring 20x6mm PTFE	2 Items	Enquire	AZLOSWN503	\$4.72
Stir Bar Pivot Ring 25x6mm PTFE	2 Items	Enquire	AZLOSWN504	\$5.07
Stir Bar Pivot Ring 35x6mm PTFE	2 Items	Enquire	AZLOSWN505	\$6.48
Stir Bar Pivot Ring 45x8mm PTFE	2 Items	Enquire	AZLOSWN506	\$7.99
Stir Bar Pivot Ring 40x8mm PTFE	2 Items	Enquire	AZLOSWN507	\$7.33
Stir Bar Pivot Ring 60x9.5mm PTFE	2 Items	Enquire	AZLOSWN508	\$12.75
Stir Bar Pivot Ring 50x8mm PTFE	2 Items	Enquire	AZLOSWN509	\$8.56
Stir Bar Pivot Ring 70x9.5mm PTFE	2 Items	Enquire	AZLOSWN510	\$14.21
Stir Bar Pivot Ring 30x6mm PTFE	2 Items	Enquire	AZLOSWN511	\$5.51
Stir Bar Oval 10x5mm PTFE	2 Items	Enquire	AZLOSWN518	\$7.19

More codes on the next page...

Bio-Strategy Ltd.

T: 0800 34 24 66

E: sales.nz@bio-strategy.com

www.bio-strategy.com

<https://shop.bio-strategy.com>

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.

Scilabware Azlon® Magnetic Stirrer and Spin Bars

Description	Pack Size	Availability	Code	Sale Price (excl. GST)
Stir Bar Oval 15x6mm PTFE	2 Items	Enquire	AZLOSWN519	\$6.63
Stir Bar Oval 20x10mm PTFE	2 Items	Enquire	AZLOSWN520	\$7.68
Stir Bar Oval 25x12mm PTFE	2 Items	Enquire	AZLOSWN522	\$10.23
Stir Bar Oval 30x16mm PTFE	2 Items	Enquire	AZLOSWN524	\$12.53
Stir Bar Oval 35x16mm PTFE	2 Items	Enquire	AZLOSWN525	\$16.33
Stir Bar Oval 40x20mm PTFE	2 Items	Enquire	AZLOSWN526	\$20.05
Stir Bar Oval 50x20mm PTFE	2 Items	Enquire	AZLOSWN528	\$31.10
Stir Bar Cylinder 6x3mm PTFE	2 Items	Enquire	AZLOSWN652	\$4.01
Stir Bar Cylinder 8x3mm PTFE	2 Items	Enquire	AZLOSWN654	\$4.07
Stir Bar Cylinder 12x4.5mm PTFE	2 Items	Enquire	AZLOSWN656	\$4.07
Stir Bar Cylinder 13x6mm PTFE	2 Items	Enquire	AZLOSWN657	\$4.47
Stir Bar Cylinder 15x4.5mm PTFE	2 Items	Enquire	AZLOSWN658	\$4.16
Stir Bar Cylinder 15x6mm PTFE	2 Items	Enquire	AZLOSWN659	\$4.54
Stir Bar Cylinder 20x6mm PTFE	2 Items	Enquire	AZLOSWN660	\$4.28
Stir Bar Cylinder 25x6mm PTFE	2 Items	Enquire	AZLOSWN662	\$4.59
Stir Bar Cylinder 25x8mm PTFE	2 Items	Enquire	AZLOSWN6620	\$6.57
Stir Bar Cylinder 30x6mm PTFE	2 Items	Enquire	AZLOSWN664	\$5.07
Stir Bar Cylinder 35x6mm PTFE	2 Items	Enquire	AZLOSWN666	\$6.09
Stir Bar Cylinder 40x8mm PTFE	2 Items	Enquire	AZLOSWN668	\$6.36
Stir Bar Cylinder 45x8mm PTFE	2 Items	Enquire	AZLOSWN670	\$7.47
Stir Bar Cylinder 50x8mm PTFE	2 Items	In stock	AZLOSWN672	\$7.81
Stir Bar Cylinder 60x9mm PTFE	2 Items	Enquire	AZLOSWN674	\$10.12
Stir Bar Cylinder 70x9.5mm PTFE	2 Items	Enquire	AZLOSWN676	\$13.19
Stir Bar Cylinder 75x13mm PTFE	2 Items	Enquire	AZLOSWN677	\$26.30
Stir Bar Cylinder 80x9.5mm PTFE	2 Items	Enquire	AZLOSWN678	\$15.18
Stir Bar Plain 108x27mm PTFE	2 Items	Enquire	AZLOSWN680	\$70.12
Stir Bar Plain 159x27mm PTFE	2 Items	Enquire	AZLOSWN682	\$87.65
Stir Bar Octagonal 25x8mm PTFE	2 Items	Enquire	AZLOSWN710	\$6.09
Stir Bar Octagonal 75x12mm PTFE	2 Items	Enquire	AZLOSWN716	\$21.67
Stir Bar Cross Head 60x17mm PTFE	2 Items	Enquire	AZLOSWN927	\$61.36
Stir Bar Cross Head 22x15mm PTFE	2 Items	Enquire	AZLOSWN933	\$18.40
Stir Bar Cross Head 30x12mm PTFE	2 Items	Enquire	AZLOSWN934	\$20.16
Stir Bar Micro 7x2mm PTFE	2 Items	Enquire	AZLOSWN950	\$5.51
Stir Bar Micro 8x1.5mm PTFE	2 Items	Enquire	AZLOSWN952	\$5.70
Stir Bar Micro 8x1.5mm PTFE	2 Items	Enquire	AZLOSWN954	\$6.40
Stir Bar Micro 13x3mm PTFE	2 Items	Enquire	AZLOSWN956	\$7.19
Stir Bar Asst Pivot Ring/Oval/Micro PTFE	Set of 22	Enquire	AZLOSWN515	\$102.91
Stirring Bar Magnetic Retriever 250mm	1 Item	Enquire	AZLOSWN564	\$19.51
Stirring Bar Magnetic Retriever 350mm	1 Item	Enquire	AZLOSWN566	\$19.36



Bio-Strategy Ltd.

T: 0800 34 24 66

E: sales.nz@bio-strategy.com

www.bio-strategy.com

<https://shop.bio-strategy.com>



bio-strategy
delivering technology

Prices & discounts above may be subject to change. Not valid with any other special offer or discounts and cannot be applied in retrospect. All offers subject to availability and while stocks last. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.