

Falcon[®] Culture Slides on Sale!

Falcon[®] Culture Slides for *in situ* analysis.

Falcon Culture Slides allow you to culture cells and then analyse them on a glass microscope slide. Cells are grown in a plastic chamber affixed to a specially prepared glass microscope slide.

Cells can be fixed and stained in place without disruption of the cell monolayer. The chamber is easily and safely removed with an easy-to-use, disposable Safety Removal Tool.

Designed for consistent cell culture results.

Speak to your local Bio-Strategy Representative for further details or call Customer Service on 0800 34 24 66 or email sales.nz@bio-strategy.com

Description	Pack Size	Code	Sale Price (excl. GST)
Falcon 8 Well Culture Slide, 12/Pack, 96/Case	96	BDAA354108	\$1,048.20
Falcon 8 Well Culture Slide, 12/Pack, 24/Case	24	BDAA354118	\$344.60
Falcon 4 Well Culture Slide, 12/Pack, 24/Case	24	BDAA354114	\$330.93
Falcon 4 Well Culture Slide, 12/Pack, 96/Case	96	BDAA354104	\$1,048.20



Starting from only \$330

**VALID UNTIL:
30 September 2019**

**Ref. Code:
BDAA-SLIDES**

Bio-Strategy Ltd.

T: 0800 34 24 66

E: sales.nz@bio-strategy.com

www.bio-strategy.com

<https://shop.bio-strategy.com>

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.