

Esco Isotherm® Laboratory Ovens

Quality range of laboratory ovens to suit your requirements.

Esco Isotherm® forced convection laboratory ovens provide controlled thermal convection drying for many laboratory applications including moisture removal and curing. Their ergonomic design, microprocessor PID control, safety protection including an independent adjustable over-temperature thermostat and high temperature alarm makes the Isotherm® oven range suitable for a multitude of laboratory applications.

- Units in or coming into stock: 54L, 110L, 170L and 240L
- Temperature range: ambient +7.5°C to 250°C
- Dual display showing set and actual temperature
- Stable heating and excellent temperature uniformity with 4-zone heated air jacket
- Fan forced air circulation with adjustable fan speed
- Superior insulation reduces external surface temperatures

Speak to your local Bio-Strategy Representative for further details or call Customer Service on 0800 34 24 66 or email sales.nz@bio-strategy.com

Description	Code	Price (excl. GST)	Sale Price (excl. GST)
Oven 54L to 250°C	ESCOOFA-54-8	\$2,730.00	\$2,185.00
Oven 110L to 250°C	ESCOOFA-110-8	\$3,570.00	\$2,856.00
Oven 170L to 250°C	ESCOOFA-170-8	\$4,725.00	\$3,780.00
Oven 240L to 250°C	ESCOOFA-240-8	\$5,890.00	\$4,712.00

Please note: all ovens feature a stainless-steel interior with two adjustable shelves.



**VALID UNTIL:
30 September 2019**

**Ref. Code:
ESCO-OVEN**

**20% Off
Current and
Incoming
Stock!**

Bio-Strategy Ltd.

T: 0800 34 24 66

E: sales.nz@bio-strategy.com

www.bio-strategy.com

<https://shop.bio-strategy.com>

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.