

Kimble Low-Form "Colorware" Beakers

VALID UNTIL:
30 September 2019

Ref. Code:
KIML-BEAKERS

"Colorware" blue beakers in stock and on sale!

Choose Kimble "Colorware" glass beakers for ease of identification in the laboratory. These beakers offer excellent mechanical strength and durability, while providing high resistance to chemical attack and thermal shock.

- Use of coloured beakers reduces the risk of cross contamination
- Easily identify glassware between labs
- All sizes feature double capacity scales with easy-to-read graduations
- Durable, colour, matte finish marking area for use with an ordinary pencil or marker
- Design meets ASTM Specification E960, Type I requirements
- Manufactured from 33 expansion, low extractable borosilicate glass conforming to USP Type I and ASTM E438, Type I, Class A requirements

Speak to your local Bio-Strategy Representative for further details or call Customer Service on 0800 34 24 66 or email sales.nz@bio-strategy.com

Description	Pack Size	Availability	Code	Sale Price (excl. GST)
Beaker, Griffin, Low-Form, Blue Scale, 1000ml	6 Items	In Stock	KIML14000B-1000	\$95.00
Beaker, Griffin, Low-Form, Blue Scale, 600ml	6 Items	In Stock	KIML14000B-600	\$48.90
Beaker, Griffin, Low-Form, Blue Scale, 250ml	12 Items	In Stock	KIML14000B-250	\$67.50



**Special
Sale Prices
Starting at
\$48.90!**

Bio-Strategy Ltd.

T: 0800 34 24 66

E: sales.nz@bio-strategy.com

www.bio-strategy.com

<https://shop.bio-strategy.com>

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.