

# Westfield Medical Autoclave Indicator Tape

**On sale, ex stock, order today!**

Autoclave tape is designed for use with flexible sterile barrier systems such as wrap and plain closure bags, although it will also adhere strongly to other surfaces such as linen. Manufactured under BS EN ISO 9001:2015 and BS EN ISO 13485:2016. Lead free.

- Printed onto a robust crepe paper-base which is coated with high tack adhesive
- Available in two standard widths of 19mm and 25mm
- Available on standard 76mm cores to allow use on tape dispenser systems
- Indicator ink on the tape conforms with ISO 11140-1:2014 for use as a Class 1 chemical indicator
- Latex free option available (please enquire)
- Sterilisation pouches & rolls available (please enquire)

**Speak to your local Bio-Strategy Representative for further details or call Customer Service on 0800 34 24 66 or email [sales.nz@bio-strategy.com](mailto:sales.nz@bio-strategy.com)**

Description	Pack Size	Code	Price (excl. GST)	Sale Price (excl. GST)
Autoclave Tape Indicating 19mmx50m	1 Roll	NZNZTAPE19MM	<b>\$8.61</b>	<b>\$6.00</b>
Autoclave Tape Indicating 25mmx50m	1 Roll	NZNZTAPE25MM	<b>\$10.89</b>	<b>\$7.50</b>



**Save  
Over  
30%**

Bio-Strategy Ltd.

T: 0800 34 24 66

E: [sales.nz@bio-strategy.com](mailto:sales.nz@bio-strategy.com)

[www.bio-strategy.com](http://www.bio-strategy.com)

<https://shop.bio-strategy.com>

**VALID UNTIL:  
30 September 2019**

**Ref. Code:  
WEST-TAPE**

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.