

OHAUS Starter 3100-F Benchtop pH Meter for Standard Laboratory Applications

Complete system including meter, electrode with ATC and flexible electrode holder.

After more than a century of perfecting the art of measurement through our durable weighing products, OHAUS precision is now available in a line of benchtop pH meters. Starter 3100 has many smart features that make it as intuitive as it is accurate for standard pH measurement. The easy-to-read large backlit LCD display informs the operator of instrument status and the auto buffer recognition feature ensures error-free operation during the calibration process.

- Applications: pH, oxidation-reduction potential (ORP) and temperature measurements
- Measurement Ranges: -2 to 16pH; -1,999 to +1,999mV; 5°C to 110°C
- Measurement Resolutions: 0.01pH; 1mV; 0.1°C
- 99 item memory and RS232 interface allows for connection to external devices
- Stand-alone electrode holder allows convenient positioning of the electrode
- 12-month warranty

Speak to your local Bio-Strategy Representative for further details or call Customer Service on 0800 34 24 66 or email sales.nz@bio-strategy.com

Description	Code	Price (excl. GST)	Sale Price (excl. GST)
Starter 3100-F pH Meter Kit	OHAU30058489	\$1,230.00	\$995.00

Includes ST3100-B meter, in-use cover, stand-alone electrode holder, ST310 pH electrode, pH buffer set (4.01, 7.01 & 10.01), power adaptor, certificate and manual.

ON SALE
Affordable pH/ORP
Measurement



Bio-Strategy Ltd.

T: 0800 34 24 66

E: sales.nz@bio-strategy.com

www.bio-strategy.com

<https://shop.bio-strategy.com>

VALID UNTIL:
30 September 2019

Ref. Code:
OHAU-3100

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.