

New PALL Universal wwPTFE Membranes

A new universal membrane from PALL.

PALL's new wwPTFE (water wettable polytetrafluoroethylene) membrane technology is ideal for aqueous solutions and offers maximum chemical compatibility for aggressive solvents. wwPTFE membrane is a low protein binding membrane which removes unwanted particulate from samples with high recovery of critical proteins.

wwPTFE Membrane Applications:

- HPLC sample preparation and clarification
- Metabolic studies
- Applications that use harsh organics
- Combinatorial chemistry library screens
- Solvent-based applications
- Low binding

Speak to your local Bio-Strategy Representative for further details or call Customer Service on 0800 34 24 66 or email sales.nz@bio-strategy.com

**VALID UNTIL:
30 September 2019**

**Ref. Code:
PALL-wwPTFE**

**10%
Off**



Description	Pack Size	Code	Sale Price (excl. GST)
13mm AcroDisc 0.2um wwPTFE Mini 300/Pk	300	PALL2400	\$965.70
13mm AcroDisc 0.2um wwPTFE Mini 1000/Pk	1000	PALL2401	\$1,796.40
13mm AcroDisc 0.45um wwPTFE Mini 300/Pk	300	PALL2402	\$965.70
13mm AcroDisc 0.45um wwPTFE Mini 1000/Pk	1000	PALL2403	\$1,796.40
OEM 13mm 0.2um wwPTFE 1000/Pk	1000	PALL6614567	\$947.70
OEM 13mm 0.45um wwPTFE 1000/Pk	1000	PALL6614563	\$947.70
25mm AcroDisc 0.2um wwPTFE 200/Cs	200	PALL4927	\$560.70
25mm AcroDisc 0.2um wwPTFE 1000/Pk	1000	PALL4929	\$1,867.50
25mm AcroDisc 0.45um wwPTFE 200/Cs	200	PALL4914	\$560.70
25mm AcroDisc 0.45um wwPTFE 1000/Pk	1000	PALL4932	\$1,867.50
ValuPrep 25mm 0.45um wwPTFE 1000/Pk	1000	PALL6604566	\$985.50
ValuPrep 25mm 0.2um wwPTFE 1000/Pk	1000	PALL6604652	\$985.50
OEM 25mm 0.2um wwPTFE 1000/Pk	1000	PALL6614562	\$824.40
OEM 25mm 0.45um wwPTFE 1000/Pk	1000	PALL6614566	\$824.40

More codes on the next page...

Bio-Strategy Ltd.

T: 0800 34 24 66

E: sales.nz@bio-strategy.com

www.bio-strategy.com

<https://shop.bio-strategy.com>

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.

New PALL Universal wwPTFE Membranes - 10% Off!

Description	Pack Size	Code	Sale Price (excl. GST)
AcroPrep Advance Filter Plate 350uL 0.2um wwPTFE 10/Pk	10	PALL8582	\$283.50
AcroPrep Advance Filter Plate 350ul 0.45um wwPTFE 10/Pk	10	PALL8584	\$283.50
AcroPrep Advance Filter Plate 1ml 0.2um wwPTFE 5/Pk	5	PALL8682	\$180.90
AcroPrep Advance Filter Plate 1ml 0.45um wwPTFE 5/Pk	5	PALL8684	\$180.90
AcroPrep Advance Filter Plate 2ml 0.2um wwPTFE 5/Pk	5	PALL8782	\$206.10
AcroPrep Advance Filter Plate 2ml 0.45um wwPTFE 5/Pk	5	PALL8784	\$206.10
wwPTFE 0.2um 25mm Disc 100/Pk	100	PALL60538	\$137.70
wwPTFE 0.2um 47mm Disc 100/Pk	100	PALL60539	\$231.30
wwPTFE 0.2um 50mm Disc 100/Pk	100	PALL60542	\$262.80
wwPTFE 0.45um 25mm Disc 100/Pk	100	PALL60547	\$137.70
wwPTFE 0.45um 47mm Disc 100/Pk	100	PALL60548	\$358.20
wwPTFE 0.45um 50mm Disc 100/Pk	100	PALL60550	\$261.00
wwPTFE 0.45um 90mm Disc 100/Pk	100	PALL60551	\$549.90
Nanosep 0.2um wwPTFE 100/Pk	100	PALLODPTFE02C34	\$321.30
Nanosep 0.2um wwPTFE 500/Pk	500	PALLODPTFE02C35	\$952.20
Nanosep 0.45um wwPTFE 100/Pk	100	PALLODPTFE04C34	\$311.40
Nanosep 0.45um wwPTFE 500/Pk	500	PALLODPTFE04C35	\$920.70
AcroDisc One PSF 0.2um wwPTFE 200/Pk	200	PALLAP-4910	\$758.70
AcroDisc One PSF 0.2um wwPTFE Autopack Tubes 200/Pk	200	PALLAP-4911	\$729.90
AcroDisc One PSF 0.2um wwPTFE 1000/Pk	1000	PALLAP-4912	\$2,648.70
AcroDisc One PSF 0.2um GxP/wwPTFE 200/Pk	200	PALLAP-4913	\$884.70
AcroDisc One PSF 0.2um GxP/wwPTFE Autopack Tubes 200/Pk	200	PALLAP-4914	\$891.00
AcroDisc One PSF 0.2um GxP/wwPTFE 1000/Pk	1000	PALLAP-4915	\$3,159.90
AcroDisc One PSF 0.45um wwPTFE 200/Pk	200	PALLAP-4916	\$758.70
AcroDisc One PSF 0.45um wwPTFE Autopack Tubes 200/Pk	200	PALLAP-4917	\$712.80
AcroDisc One PSF 0.45um wwPTFE 1000/Pk	1000	PALLAP-4918	\$2,528.10
AcroDisc One PSF 0.45um GxP/wwPTFE 200/Pk	200	PALLAP-4919	\$884.70
AcroDisc One PSF 0.45um GxP/wwPTFE Autopack Tubes 200/Pk	200	PALLAP-4920	\$884.70
AcroDisc One PSF 0.45um GxP/wwPTFE 1000/Pk	1000	PALLAP-4921	\$3,159.90



Bio-Strategy Ltd.

T: 0800 34 24 66

E: sales.nz@bio-strategy.com

www.bio-strategy.com

<https://shop.bio-strategy.com>



bio-strategy
delivering technology

Prices & discounts above may be subject to change. Not valid with any other special offer or discounts and cannot be applied in retrospect. All offers subject to availability and while stocks last. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.