

NGS Sample Quality Control Using the QIAxcel Advanced System

Start-up package includes all your QIAxcel reagents for NGS library QC analysis for up to 1,200 samples.

The revolutionary QIAxcel Advanced system replaces traditional, labour-intensive gel analysis of DNA and RNA, streamlining workflows and reducing time to result.

- Rapid analysis of up to 96 samples without manual intervention
- Rapid analysis of 12 samples in 13 minutes
- Safety and convenience with ready-to-use gel cartridges
- Accurate analysis of low concentration nucleic acids with 3–5bp resolution
- User-friendly analysis software that supports electronic signatures and records
- Objective RNA quality measurement with new RNA Integrity Score (RIS)

Speak to your local Bio-Strategy Representative for further details or call Customer Service on 0800 34 24 66 or email sales.nz@bio-strategy.com

Description	Code
QIAxcel Advanced System*	QIAG9001941
QIAxcel DNA High Resolution Kit (1200)	QIAG929002
QIAxcel Alignment Marker 15 bp/3 kb (1.5ml)	QIAG929522
QIAxcel Size Marker 50-800 bp (50ul)	QIAG929561

**Only NGS related kits mentioned. Please enquire about other kits available for DNA and RNA analysis.*

**VALID UNTIL:
30 September 2019**

**Ref. Code:
QIAxcel**

**Special
Package Price
\$48,300**



Bio-Strategy Ltd.

T: 0800 34 24 66

E: sales.nz@bio-strategy.com

www.bio-strategy.com

<https://shop.bio-strategy.com>

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.