

PHCbi Large Capacity Upright -86°C Freezers

VALID UNTIL:
30 September 2019

Ref. Code:
CASCADE

Choose between 728-litre or 845-litre models to suit your requirements.

The PHCbi (formerly Panasonic/Sanyo) MDF-U74V and MDF-DU900V-PE continue the legacy of Ultra-Low Temperature -86°C freezers with vacuum panel insulation technology (VIP) and Cool Safe Compressors to provide reliable units with increased compressor longevity and lower heat output. The units are quiet in operation, run off a standard 10-amp supply and the easy access airflow filter reduces maintenance times.

- Cascade refrigeration with Cool Safe compressors – proven reliable cooling technology with compressors running at lower temperatures
- VIP Plus insulation – provides up to 30% more storage capacity compared to conventionally insulated freezers with the same footprint
- User friendly eye-level LED display with audible and visual alarms
- Two insulated inner doors with eagle-latches reduces heat loss
- Quiet in operation – 51dB(A) and low power consumption
- Keyed door lock for security

Speak to your local Bio-Strategy Representative for further details or call Customer Service on 0800 34 24 66 or email sales.nz@bio-strategy.com

Heavy Price Reduction Applies to Stocked Units Only!

Description	Code	Price (excl. GST)	Sale Price (excl. GST)
Freezer, Upright, -86°C, 728L	SANYMDF-U74V	\$20,841.00	\$12,995.00
Freezer, Upright, -86°C, 845L	SANYMDF-DU900VX-PE	\$25,879.00	\$16,350.00

Supplied with three adjustable stainless-steel shelves.



Bio-Strategy Ltd.

T: 0800 34 24 66
E: sales.nz@bio-strategy.com
www.bio-strategy.com
<https://shop.bio-strategy.com>

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.