

## Scharlau Aquagent® Pyridine-Free Karl Fischer Reagents for Volumetric Water Determination

### Safer, faster, more stable and on sale!

Aquagent® reagents contain non-toxic Imidazole instead of the toxic Pyridine as its base. Imidazole has a better buffer capacity, giving you a faster reaction, stable end-point and more reproducible results.

**One Component System** - Simple to use but must be re-titrated frequently due to the reactivity of their components.

**Two Component System** - The solvent-component doesn't just act as a solvent medium but contains part of the reagents for a longer shelf-life and doesn't need re-titration.

**Speak to your local Bio-Strategy Representative for further details or call Customer Service on 0800 34 24 66 or email [sales.nz@bio-strategy.com](mailto:sales.nz@bio-strategy.com)**

Our most commonly supplied reagents are as follows, if you need something else please ask

Description	Pack Size	Availability	Code	Sale Price (excl. GST)
<b>Aquagent® Complet 5 KF Reagent</b> General purpose reagent for samples with high and medium water content. It has a titre of approx. 5mg water/ml. It is generally used in conjunction with methanol as a solvent.	1 Litre	In Stock	SCARAQ00031000	<b>\$81.60</b>
<b>Aquagent® Titrant 2 Pyridine Free</b> General purpose reagent that contains iodine and methanol. Titre is approx. 5mg water/ml. Must be used in conjunction with Aquagent® Solvent.	1 Litre	In Stock	SCARAQ00061000	<b>\$88.90</b>
<b>Aquagent® Titrant 5 Pyridine Free</b> General purpose reagent that contains iodine and methanol. Titre is approx. 2mg water/ml. Must be used in conjunction with Aquagent® Solvent.	1 Litre	In Stock	SCARAQ00011000	<b>\$69.20</b>
<b>Aquagent® Solvent Pyridine Free</b> General purpose reagent that contains SO <sub>2</sub> , imidazole and methanol. It must be used in conjunction with Aquagent® Titrant.	1 Litre	In Stock	SCARAQ00021000	<b>\$87.80</b>
<b>Aquagent® Solvent CM Pyridine Free</b> Solvent-component for titration of fats and oils. It contains chloroform, which improves solubility of long-chained hydrocarbons.	1 Litre	In Stock	SCARAQ00081000	<b>\$93.90</b>

**STOCK UP  
NOW AND  
SAVE**

**VALID UNTIL:  
30 September 2019**

**Ref. Code:  
AQUAGENT**



Bio-Strategy Ltd.

T: 0800 34 24 66

E: [sales.nz@bio-strategy.com](mailto:sales.nz@bio-strategy.com)

[www.bio-strategy.com](http://www.bio-strategy.com)

<https://shop.bio-strategy.com>

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.