

# VWR Everyday Consumables

## Economical solutions for your laboratory.

From autoclave bags to PCR plasticware – Bio-Strategy exclusively offers a large range of quality consumables from VWR.

Speak to your local Bio-Strategy Representative for further details or call Customer Service on 0800 34 24 66 or email [sales.nz@bio-strategy.com](mailto:sales.nz@bio-strategy.com)

### Sterile spreaders, loops and spatulas

Economical sterile tools to help in microbiology, food analysis and cell culture.

Description	Pack Size	Availability	Code	Sale Price (excl. GST)
Spreaders, L-Shaped, Disposable, Blue, Sterile 5/Bag	1000	Stocked	VWRI612-1561	<b>\$154.50</b>
Spreaders, L-Shaped, Disposable, Blue, Sterile Individually Packed	500	Enquire	VWRI612-1560	<b>\$181.51</b>
Inoculating Loops Hard 10µl, Disposable, Blue, Sterile 20/Bag	1000	Enquire	VWRI612-9354	<b>\$82.80</b>
Inoculating Loops Hard 1µl, Disposable, Green, Sterile 20/Bag	1000	Enquire	VWRI612-9352	<b>\$82.80</b>
Inoculating Loops Soft 1µl, Disposable, Light Green, Sterile 20/Bag	500	Enquire	VWRI612-9355	<b>\$42.70</b>
Inoculating Loops Soft 10µl, Disposable, Light Blue, Sterile 10/Bag	500	Enquire	VWRI612-9357	<b>\$42.70</b>
Inoculating Loops Soft 10µl, Disposable, Light Blue, Sterile 20/Bag	1000	Enquire	VWRI612-9358	<b>\$82.80</b>
Spatulas Disposable Opaque, 210mm/7mm	300	Enquire	VWRI231-0104	<b>\$69.23</b>
Spatulas Disposable Natural, 210mm/7mm, Sterile and Non-Pyrogenic	100	Enquire	VWRI231-0593	<b>\$295.05</b>

### Disposable cuvettes (cells)

Disposable cuvettes for 300-900nm spectrophotometers semi-micro and macro, PS or PMMA.

Description	Pack Size	Availability	Code	Sale Price (excl. GST)
Cuvette Macro, PS, 340-900nm, Volume 2.5-4.5ml, 10mm	100	Stocked	VWRI634-0675	<b>\$12.64</b>
Cuvette Macro, PMMA, 300-900nm, Volume 2.5-4.5ml, 10mm	100	Stocked	VWRI634-0677	<b>\$20.11</b>
Cuvette Semi-Micro, PS, 340-900nm, Volume 1.5-3ml, 10mm	100	Stocked	VWRI634-0676	<b>\$13.49</b>
Cuvette Semi-Micro, PMMA, 300-900nm, Volume 1.5-3ml, 10mm	100	Enquire	VWRI634-0678	<b>\$20.18</b>

**More products on promotion on the next page...**

Bio-Strategy Ltd.

T: 0800 34 24 66

E: [sales.nz@bio-strategy.com](mailto:sales.nz@bio-strategy.com)



[www.bio-strategy.com](http://www.bio-strategy.com)

<https://shop.bio-strategy.com>

**VALID UNTIL:**  
**30 September 2019**

**Ref. Code:**  
**VWR-LAB**

**Save  
over  
30%**



Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.

## Filter paper, glass microfibre filters and extraction thimbles

ON SALE

### Quantitative filter paper:

- 454: Retention of zirconium hydroxide, aluminium hydroxide, cobalt sulphide and other coarse or gelatinous precipitates
- 474: General quantitative analyses, also for filtration of aluminium phosphate, zinc ammonium phosphate and silver chromate

### Qualitative filter paper:

- 415 and 301: Coarse precipitates such as aluminium hydroxide, cobalt carbonate, mercury iodide, magnesium hydroxide and more
- 417: Rapid filtrations where fine particle retention is not required
- 600 and 310: Qualitative filter paper for use in routine laboratories, universities and hospitals, made of 100% cellulose with no binder

### Glass fibre filters:

- Grade 691: Ideal for radioimmunoassay of weak beta emitters by scintillation counting and for filtration of precipitated proteins, cells and cell debris, also used for gravimetric determination of airborne particulates and air pollution monitoring
- Grade 693: Used for clarification and monitoring of suspended solids in wastewater and water and in marine chemistry to filter particulate components in sea water and in biochemistry for cell harvesting in radioimmunoassay procedures



Description	Pack Size	Availability	Code	Sale Price (excl. GST)
Filter Paper Qualitative 415, Fast, 12-15µm Particle Retention, 185mm	100	Stocked	VWRI516-0836	\$20.98
Filter Paper Qualitative 415, Fast, 12-15µm Particle Retention, 250mm	100	Stocked	VWRI516-0838	\$32.88
Filter Paper Quantitative 454, Fast, 12-15µm Particle Retention, 185mm	100	Stocked	VWRI516-0857	\$58.62
Filter Paper Qualitative 417, Rapid, 40µm Particle Retention, 90mm	50	Enquire	VWRI516-0821	\$8.45
Filter Paper Qualitative 417, Rapid, 40µm Particle Retention, 110mm	50	Enquire	VWRI516-0822	\$8.96
Filter Paper Qualitative 417, Rapid, 40µm Particle Retention, 125mm	50	Enquire	VWRI516-0823	\$10.25
Filter Paper Quantitative 474, Medium, 5-10µm Particle Retention, 90mm	100	Enquire	VWRI516-0847	\$21.42
Filter Paper Quantitative 474, Medium, 5-10µm Particle Retention, 110mm	100	Enquire	VWRI516-0848	\$27.44
Filter Paper Qualitative 413, Medium, 5-13µm Particle Retention, 150mm	100	Enquire	VWRI516-0817	\$12.54
Filter Paper Qualitative 600 (No Binder) 420x520mm	500	Enquire	VWRI516-0309	\$107.18
Glass Fibre Filters 693, 1.2µm Particle Retention, 24mm	100	Enquire	VWRI516-0868	\$20.40
Glass Fibre Filters 693, 1.2µm Particle Retention, 47mm	100	Enquire	VWRI516-0869	\$29.55
Glass Fibre Filters 691, 1.6µm Particle Retention, 150mm	100	Stocked	VWRI516-0867	\$124.85
Glass Fibre Filters 691, 1.6µm Particle Retention, 125mm	100	Stocked	VWRI516-0866	\$103.94
Cellulose Extraction Thimbles 501, 25x80mm	25	Stocked	VWRI516-0253	\$91.76
Cellulose Extraction Thimbles 501, 33x80mm	25	Enquire	VWRI516-0259	\$105.50
Cellulose Extraction Thimbles 501, 22x80mm	25	Enquire	VWRI516-0252	\$99.50

## PCR tubes, strips and plates

A vast range of performance tubes for various centrifugation speeds.

Description	Pack Size	Availability	Code	Sale Price (excl. GST)
8-Tube Strips for PCR, with Individually Attached Domed Caps, Clear, 0.2ml	120	Stocked	VWRI732-0545	<b>\$103.80</b>
Individual PCR Tubes, with Attached Flat Caps, 0.2ml	1000	Enquire	VWRI732-0548	<b>\$46.50</b>
8-Tube Strips, Low Profile, Qpcr, Ind. Attached Flat Caps, Clear, 0.2ml	120	Stocked	VWRI211-0339	<b>\$137.90</b>
8-Tube Strips for PCR, with Attached Domed Cap Strips, 0.2ml	125	Stocked	VWRI732-0546	<b>\$139.80</b>
PCR Plates 96-Well Standard Profile, Non-Skirted, Alphanumeric, 0.2ml	100	Stocked	VWRI211-0262	<b>\$287.50</b>
PCR Plates 96-Well Standard Profile, Non-Skirted, Alphanumeric, 0.2ml, PP	100	Enquire	VWRI732-2387	<b>\$319.80</b>



ON SALE

## Autoclave bags

Quality autoclave and biohazard bags for 135-145°C sterilisation. Orange models contain a sterilisation indicator. Sturdy blender bags for use with homogenisers.

Description	Pack Size	Availability	Code	Sale Price (excl. GST)
Bag Autoclavable Clear, Print Biohazard, PP, 145°C, 305x660mmx40µm	200	Stocked	VWRI129-0587	<b>\$37.06</b>
Bag Autoclavable Clear, Print Biohazard, PP, 145°C, 600x780mmx40µm	200	Stocked	VWRI129-0589	<b>\$89.39</b>
Bag Autoclavable, Plain, PP, 145°C, 900x1100mmx40µm	100	Stocked	VWRI129-0586	<b>\$92.62</b>
Bag Autoclavable, Plain, PP, 145°C, 700x1100mmx40µm	350	Enquire	VWRI129-0585	<b>\$205.80</b>
Bag Autoclavable, Plain, PP, 145°C, 420x610mmx40µm	500	Enquire	VWRI129-0582	<b>\$129.20</b>
Bag Autoclavable, Clear Print Biohazard, PP, 145°C, 420x610mmx40µm	200	Stocked	VWRI129-0588	<b>\$45.69</b>
Bag Autoclavable, Print Biohazard Orange w. Indicator PP, 483x610mmx50µm	200	Enquire	VWRI129-0025	<b>\$244.50</b>
Bag Autoclavable, Print Biohazard Orange w. Indicator PP, 635x889mmx50µm	200	Enquire	VWRI129-0026	<b>\$443.50</b>



ON SALE

Bio-Strategy Ltd.

T: 0800 34 24 66

E: [sales.nz@bio-strategy.com](mailto:sales.nz@bio-strategy.com)

[www.bio-strategy.com](http://www.bio-strategy.com)

<https://shop.bio-strategy.com>



**bio-strategy**  
delivering technology

Prices & discounts above may be subject to change. Not valid with any other special offer or discounts and cannot be applied in retrospect. All offers subject to availability and while stocks last. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.