

# Axygen New Long Reach Pipette Tips

**VALID UNTIL:**  
31 March 2020

**Ref. Code:**  
AXYG-TIPS

## Special introductory price - ideal for pipetting from long and narrow tubes


Axygen's NEW Long Reach Tips are longer and thinner to help you minimise contamination when aspirating from long and narrow storage tubes, cell culture flasks, plates and sample bags. The tips are elongated, clear and compatible with the vast majority of research grade pipettors.

- Slender Elongated Tip
- Graduation marks for easy visualisation of liquid volume
- DNase, RNase, pyrogen free
- Autoclavable (121°C for 15 minutes) or purchase already sterilised tips
- 10µl tip length up to 45mm
- 200µl tip length up to 77.5mm
- 1000µl tip extended length (102mm)

**From  
20%  
OFF**

**Speak to your local Bio-Strategy Representative for further details or call Customer Service on 0800 34 24 66 or email [sales.nz@bio-strategy.com](mailto:sales.nz@bio-strategy.com)**

Description	Unit of Measure	Code	List Price (excl. GST)	Sale Price (excl. GST)
Axygen 10µl long reach tips, bulk, non-sterile	1 * 1000 items	AXYGT-10-XT-C	<b>\$45.77</b>	<b>\$23.58</b>
Axygen 10µl long reach tips, rack, sterile	1 * 960 items	AXYGT-10XT-C-R-S	<b>\$67.54</b>	<b>\$50.78</b>
Axygen 200µl long reach tips, bulk, non-sterile	1 * 1000 items	AXYGT-200XT-C	<b>\$48.92</b>	<b>\$38.44</b>
Axygen 200µl long reach tips, rack, sterile	1 * 960 items	AXYGT-200XT-C-R-S	<b>\$82.06</b>	<b>\$51.40</b>
Axygen 1000µl long reach tips, bulk, non-sterile	1 * 768 items	AXYGT-1000XT-C	<b>\$44.82</b>	<b>\$23.08</b>
Axygen 1000µl long reach tips, rack, sterile	1 * 768 items	AXYGT-1000XT-C-R-S	<b>\$73.86</b>	<b>\$42.80</b>



Bio-Strategy Ltd.

T: 0800 34 24 66

E: [sales.nz@bio-strategy.com](mailto:sales.nz@bio-strategy.com)

[www.bio-strategy.com](http://www.bio-strategy.com)

<https://shop.bio-strategy.com>

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.