

Corning® 15mL & 50mL PP Centrifuge Tubes

VALID UNTIL:
31 March 2020

Ref. Code:
CORNCENT-PROMO

Corning® 15ml & 50ml Clear Polypropylene (PP) Centrifuge Tubes Come with The Advanced Polyethylene Centristar™ Caps and Are Bulk Packed. Foam Racks Can Be Purchased Separately. The Centristar™ Cap Has an Easy-On/Easy-Off Flat Top and Offers Advanced Ergonomics with Its Wider Knurls And Roll-Over Edge Design For Easier Gripping. This Leak-Proof Design Comes with A Revolutionary Plug Feature That Virtually Eliminates All Seepage When Used Under Recommended Conditions.

Features include:

- Maximum RCF is 12,500 xg
- Temperature Range From -80°C To 120°C
- 5 Year Shelf Life

Speak to your local Bio-Strategy Representative for further details or call Customer Service on 0800 34 24 66 or email sales.nz@bio-strategy.com

**Corning
Centrifuge
Tube Sale**

Ref. Code: CORNCENT-PROMO

Valid until 31/03/2020

Description	Pack Size	Code	Sale Price (excl. GST)
Corning® 15mL PP Centrifuge Tubes, Bulk Packed with CentriStar™ Cap, Sterile, 500/Case	500	CORN430791	\$93.00
Corning® 50mL PP Centrifuge Tubes, Conical Bottom with CentriStar™ Cap, Bulk Packed, Sterile, 25/Sleeve, 500/Case	500	CORN430829	\$117.00



Bio-Strategy Ltd.

T: 0800 34 24 66
E: sales.nz@bio-strategy.com
www.bio-strategy.com
<https://shop.bio-strategy.com>

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.