

# Esco Airstream E 1.2 Class II Microbiological Safety Cabinet

**Protection you can count on: operator, product and environment protection!**

Esco offers an excellent range of Class II Cabinets, the Airstream E range continues this tradition. The Esco Airstream E Class II Microbiological Safety Cabinets are available in widths from 0.6 metres to 1.8 metres. Tempered glass side panels allow entry of natural light and the tempered glass angled manual sash provides an uninhibited view of the work zone. The 1.2m Airstream E unit offered includes 2xpower outlets, fluorescent and UV lighting as standard. A full width armrest ensures correct posture with unimpeded air inflow to the cabinet.

- Quiet: only 53.5dB(A) and energy efficient: only 160W
- ULPA downflow and exhaust filters with 99.999% efficiency, 10X better than HEPA filter (99.99%) improves safety and delivers ISO Class 3 air cleanliness
- Built-in UV lamp with timer safety interlocked to prevent UV exposure
- Sentinel Gold microprocessor controller with digital display shows cabinet status
- Multi-piece stainless-steel work surface is easy to access and clean
- Night Setback Mode reduces energy consumption by a further 60%

**Speak to your local Bio-Strategy Representative for further details or call Customer Service on 0800 34 24 66 or email [sales.nz@bio-strategy.com](mailto:sales.nz@bio-strategy.com)**

**VALID UNTIL:  
31 March 2020**

**Ref. Code:  
ESCO-AIR**

**Safety at an  
Affordable  
Price!**

Ref. Code: ESCO-AIR

Valid until 31/03/2020

Description	Code	Price (excl. GST)	Sale Price (excl. GST)
Esco Airstream E 1.2m Class II Cabinet	ESCOAC2-4E8	<b>\$13,793.00</b>	<b>\$9,975.00</b>

*Please note: sales price includes 1.2m Class II Cabinet with 2xelectrical outlets, 1xUV and a free adjustable height stand with lockable castors (valued at \$1,172.00).*



Bio-Strategy Ltd.

T: 0800 34 24 66  
E: [sales.nz@bio-strategy.com](mailto:sales.nz@bio-strategy.com)  
[www.bio-strategy.com](http://www.bio-strategy.com)  
<https://shop.bio-strategy.com>

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.