

## Gyrozen GZ-0624 Centrifuge

### The Smallest Refrigerated Centrifuge to Spin 15, 50 or 100ml Samples

- For up to 6,000 Rpm (3,960 Xg) With 6 x 50 ml Conical Tubes and 4,800 Rpm (3,888 xg) With 4 x 100 ml Tubes
- Commonly used in the Fractionation Of Blood Samples and Simple Precipitation of Biological Samples
- Fast Cooling Function to 4°C in 5 Minutes For Fast Startup Of Cooled Samples
- Three Kinds of Fixed Angle Rotors for Micro Tubes, 15 ml And 50 ml Tubes As Well As Four Swing-Out Rotors.

Speak to your local Bio-Strategy Representative for further details or call Customer Service on 0800 34 24 66 or email [sales.nz@bio-strategy.com](mailto:sales.nz@bio-strategy.com)

**VALID UNTIL:**  
31 March 2020  
**Ref. Code:**  
GYROGZ-0624R

**10% OFF**

Description	Code	UoM	Price (excl. GST)	Sale Price (excl. GST)
624R, Low-Speed Refrigerated Centrifuge without Rotor (220V, 50/60Hz)	GZ-0624R-B	1	<b>\$8,361.00</b>	<b>\$6,600.00</b>
Four-wing Swing 50ml rotor of 624R	GRS-SR-50-4	1	<b>\$612.00</b>	<b>\$551.00</b>
Conical 50 ml Bucket to accommodate 1 x 50 ml conical tube	GLB-C50-50	4	<b>\$124.00</b>	<b>\$117.00</b>



Purchase from  
our webshop!

Bio-Strategy Ltd.

T: 0800 34 24 66  
E: [sales.nz@bio-strategy.com](mailto:sales.nz@bio-strategy.com)  
[www.bio-strategy.com](http://www.bio-strategy.com)  
<https://shop.bio-strategy.com>

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.