

# Pall MicroCheck® II Beverage Monitor for Microbial Analysis!

## Easy-to-use Disposable Filter Funnels to Meet Microbial Analysis Needs

- Can convert to a Petri dish for culturing or easily remove the membrane for placing on separate agar dish.
- Ready to use; simply remove from box and place directly on a filter manifold. No need for rubber stoppers or adapters.
- Volume graduations marked clearly for easy reading from any point of view.
- Volume graduations marked externally on the funnel to prevent potential hold-up of organisms.

#### APPLICATION:

- Finished product and raw material testing for microbial contamination

**Speak to your local Bio-Strategy Representative for further details or call Customer Service on 1800 008 453 or email [sales.au@bio-strategy.com](mailto:sales.au@bio-strategy.com)**

**VALID UNTIL:  
31 March 2020**

**Ref. Code:  
PALLFUNNEL-PROMO**

**On Sale!**

**Ref. Code: PALLFUNNEL-PROMO**

**Valid until 31/03/2020**

Description	Pack Size	Code	Sale Price (excl. GST)
Microcheck 100mL 0.45um GN6 Metrical Membrane Gridded	50	PALL4761	<b>\$130.00</b>



Bio-Strategy Ltd.

T: 0800 34 24 66

E: [sales.nz@bio-strategy.com](mailto:sales.nz@bio-strategy.com)

[www.bio-strategy.com](http://www.bio-strategy.com)

<https://shop.bio-strategy.com>

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.