

# PHCbi MLR-352H Upright Plant Growth Chamber, 294 Litres On Sale

**VALID UNTIL:**  
31 March 2020

**Ref. Code:**  
MLR-352H

**Versatile upright plant growth chamber with flexible programmability of temperature, lighting and humidity.**

- Temperature range: 0 to 50°C with lights off. +5 to 50°C with lights on
- Lighting: 0 lux to 20,000 lux via 15x40-watt fluorescent tubes in 6 increments: 0, 1, 2, 3, 4 & 5
- Humidity: 60%RH to 90%RH (+15 to +45°C) lights off, 55%RH to 85%RH (+15 to +45°C) lights on

The unit has storage for 10 programs. Each program can have up to 12-steps with repeat 1-98 times or unlimited. Alarms for temperature (upper and lower), humidity (upper and lower), sensor, filter and door open.

- Graphic LCD for easy operation
- Programmable temperature, lighting and humidity functions
- Microprocessor PID control and refrigeration capacity control
- Data logging function
- Automatic (3-timer pattern) and manual defrost
- 4 adjustable and 1 fixed (bottom) shelves


**Speak to your local Bio-Strategy Representative for further details or call Customer Service on 0800 34 24 66 or email [sales.nz@bio-strategy.com](mailto:sales.nz@bio-strategy.com)**

**Flexible programming of temperature, and more.**

**Ref. Code: MLR-352H**

**Valid until 31/03/2020**

Description	Code	List Price (excl. GST)	Sale Price (excl. GST)
Upright Plant Growth Chamber, 294 litres	SANYMLR-352H-PE	<b>\$22,905.00</b>	<b>\$18,975.00</b>



Bio-Strategy Ltd.

T: 0800 34 24 66

E: [sales.nz@bio-strategy.com](mailto:sales.nz@bio-strategy.com)

[www.bio-strategy.com](http://www.bio-strategy.com)

<https://shop.bio-strategy.com>

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.