

VWR Filter Paper and Glass Fibre Filters on Special!

VALID UNTIL:
31 March 2020

VWR Filter paper and Glass Fibre Filters for every application!

Ref. Code:
VWR-FILTER

Qualitative:

Grade 413 and 303: For routine laboratory work, suitable for student requirements in general chemistry, also for filtration of aluminium phosphate, barium carbonate, lead chromate and more.

Grade 415 and 301: Coarse precipitates such as aluminium hydroxide, cobalt carbonate, mercury iodide, magnesium hydroxide and more.

Grade 417: Rapid filtrations where fine particle retention is not required.

Quantitative:

454: Retention of zirconium hydroxide, aluminium hydroxide, cobalt sulphide and other coarse or gelatinous precipitates.

474: General quantitative analyses and for filtration of aluminium phosphate, zinc ammonium phosphate and silver chromate.

Glass fibre filters:

Grade 691: Ideal for radioimmunoassay of weak beta emitters by scintillation counting and for filtration of precipitated proteins, cells and cell debris. Also used for gravimetric determination of airborne particulates and air pollution monitoring.

Grade 693: Used for clarification and monitoring of suspended solids in wastewater and water and in marine chemistry to filter particulate components in sea water and in biochemistry for cell harvesting in radioimmunoassay procedures and for binders analysis.

More grades and sizes available on demand!

Speak to your local Bio-Strategy Representative for further details or call Customer Service on 0800 34 24 66 or email sales.nz@bio-strategy.com

Economic Filter Solutions!

Description	Pack Size	Code	Sale Price (excl. GST)
Extraction thimbles 501 25 X 80 MM	[25 ITEMS]	VWRI516-0253	\$104.34
Filterpaper Qualitative 413 90MM	[100 ITEMS]	VWRI516-0814	\$8.11
Filterpaper Qualitative 417 90MM	[50 ITEMS]	VWRI516-0821	\$8.40
Filterpaper Qualitative 417 110MM	[50 ITEMS]	VWRI516-0822	\$9.25
Filterpaper Qualitative 415 185MM	[100 ITEMS]	VWRI516-0836	\$23.87
Filterpaper Qualitative 415 250MM	[100 ITEMS]	VWRI516-0838	\$37.39
Filterpaper Quantitative 474 90MM	[100 ITEMS]	VWRI516-0847	\$22.11
Filterpaper Quantitative 474 110MM	[100 ITEMS]	VWRI516-0848	\$28.33
Filterpaper Quantitative 454 185MM	[100 ITEMS]	VWRI516-0857	\$66.67
Glass microfibre filter 691 125MM	[100 ITEMS]	VWRI516-0866	\$118.20
Glass microfibre filter 691 150MM	[100 ITEMS]	VWRI516-0867	\$141.97
Glass microfibre filter 693 55MM	[100 ITEMS]	VWRI516-0870	\$34.83
Glass microfibre filter 693 Ø 110MM	[100 ITEMS]	VWRI516-0888	\$85.27


Bio-Strategy Ltd.

T: 0800 34 24 66

E: sales.nz@bio-strategy.com

www.bio-strategy.com

<https://shop.bio-strategy.com>



Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.