

# Leica DM750 Microscope

VALID UNTIL:  
30 September 2020

Ref. Code:  
LEIM-DM750

## For high quality images and reliable performance!

**Leica Optics + Bright LED illumination leads to a higher quality image.**

Based on the same optical platform as Leicas' research microscopes.

**Quality materials & construction gives reliable performance over a longer life.**

Weighted focus knobs with internal brass mechanism for accurate focus capability.

Wear resistant, special stage finish offers additional protection from friction damage.

AGTREAT™ on the microscope touchpoints inhibits bacteria growth and helps prevent the spread of disease from person to person.

**Add the ICC50W  
Camera module!**



### Would sharing, capturing and archiving your microscope images be useful?

Add: ICC50W High Definition Wireless Camera module

5 modes/choices of connectivity and working with images:

1. WiFi mode + Leica AirLab App for mobile devices
2. USB mode connect to PC & Leica software
3. HDMI port direct to HD screens
4. SD card
5. Ethernet

Description	Code	Price	Sale Price
Leica DM750 microscope basic outfit. Brightfield only. 45°EZ tube with fixed eyepieces, Plan 4x, 10x, 40x and 100x oil objectives	13613403	\$2,410	
Leica ICC50W camera Module	13613735	\$2,533	
Leica DM750 Microscope Kit: Brightfield & Phase contrast (Phase slider). Koehler illumination. 5 hole nosepiece. 30° binocular tube with focusing eyepieces, HI Plan 4x, 10x, 20x, 40x, & N plan 100x oil Phase objectives	LEIM-DM750	\$8,412	\$5,468
DM750 customised for your needs	Customised	POA	

Bio-Strategy Ltd.

T: 0800 34 24 66  
sales.nz@bio-strategy.com  
www.bio-strategy.com  
shop.bio-strategy.com