


New QIAcuity PCR System from QIAGEN

The QIAcuity Digital PCR System makes absolute quantification accessible and affordable for all laboratories. The walk-away automation integrates and streamlines the entire digital PCR workflow of partitioning, thermocycling and imaging into a single instrument with minimal hands-on time. It is also easy to adapt your current qPCR assays to the QIAcuity Digital PCR System. No change in plate handling is required when coming from qPCR, assuring fast assay setup and quick results.

The QIAcuity is designed to deliver precise and multiplexed quantification results for mutation detection, copy number variation (CNV), gene expression studies, gene-editing analysis, and many more.

- Fully integrated system
- Flexible sample throughput
- Comprehensive results in under 2 hours
- Scalable format (1-, 4- and 8-plate instruments)
- Advanced multiplexing capabilities (up to 5 plex)

**Demos will
be available
soon!**



Bio-Strategy Ltd.