

Esco Airstream E Class II Microbiological Safety Cabinet

VALID UNTIL:
31 December 2020

Ref. Code:
ESCO-AIR

Offering operator, environment, and product protection – just what is needed in today’s market.

The Esco AC2-4E8 Airstream E 1.2m Class II Microbiological Safety Cabinet includes a safety interlocked UV lamp with timer for secondary decontamination, 2 electrical sockets, glass side panels allowing entry of natural light (excellent for viewing notes while working in the chamber!), multi-piece work surface - easily removable for cleaning/autoclaving and an adjustable height stand with lockable castors. Operation is easy – just raise ergonomic sloping sash to operating height marks, lights and fan will activate. After 3 minutes cabinet will advise of stable airflow = ready for use.

**25%
OFF!**

- Real-time airflow sensor with alarms ensures your safety
- One of the most energy efficient cabinets available today – 70% energy savings when compared to Class II cabinets with AC motors
- Programmable safety interlocked UV light with timer – extends UV lamp life
- ULPA downflow and exhaust filters – 10 x filtration efficiency of HEPA filters creating ISO Class 3 work zone (instead of industry standard ISO Class 5)
- Quiet – only 53.5dB(A)
- Selectable Quickstart mode for fast operation and Night Setback mode for reduced power consumption by up to 60%



Description	Code	Price (ex GST)	Sale Price (ex GST)
Airstream E 1.2m Class II Type A2 Microbiological Safety Cabinet with UV Light, Electrical Outlets and Stand	ESCOAC2-4E8	\$14,822.00	\$10,755.00

Bio-Strategy Ltd.

T: 0800 34 24 66
sales.nz@bio-strategy.com
www.bio-strategy.com
shop.bio-strategy.com