

Esco Ascent™ Opti Ductless Fume Hood

VALID UNTIL:
31 December 2020

Ref. Code:
ESCO-FUME

Package deal available today from Bio-Strategy.

The Esco Ascent™ Opti Ductless Fume Hood provides both operator and environment protection by adsorbing chemical fumes onto the Nanocarb™ activated carbon filter. We are offering the ADC-3A1 hood including the optional CFA-6 Type A Nanocarb™ filter as this filter is designed for all common chemicals (toluene, xylene and many other organic compounds) found in laboratories including pathology laboratories. Options include Nanocarb filters for trapping other compounds and mobile cart.

- Angled hinged acrylic front panel with hand slots allowing easy access to work zone
- Acrylic side panels and front panel for superior viewing
- Spill retaining recessed stainless-steel work surface
- Integrated baffle improves containment and ensures the efficient removal of fumes from the work zone
- Advanced Sentinel™ Silver microprocessor controller with digital display and airflow alarm
- Fluorescent lighting and pass-through ports for electrical or service connections

**10%
Package
Discount!**



Description	Code	Price (ex GST)	Sale Price (ex GST)
Ascent Opti Ductless Fume Hood, 0.9m	ESCOSPD-3A1	\$3,892.00	
Nanocarb Standard Filter Type A CFA-6	ESCO5080013	\$953.00	
		(Total) \$4,845.00	\$4,360.00

Bio-Strategy Ltd.

T: 0800 34 24 66
sales.nz@bio-strategy.com
www.bio-strategy.com
shop.bio-strategy.com