

# Grant Stirred Waterbath Model T100-ST26

VALID UNTIL:  
31 December 2020

Ref. Code:  
GRAN-BATH

## 26-litre bath with immersion thermostat.

The Grant T100 series of general-purpose heating circulating baths offer accurate temperate control for straightforward applications in the range of 0°C to 100°C. The T100 heating circulator can be used in combination with any Grant stainless steel or plastic tank or with a T-CLAMP for independent use of virtually any vertical sided tank. Our offer has the T100 combined to a 26-litre stainless steel tank with bridge piece.

- Heating range ambient +5°C to 100°C (below ambient work possible with cooling unit)
- Stability  $\pm 0.05^\circ\text{C}$
- Three programmable temperature pre-sets
- Low liquid protection and fixed over temperature cut-out
- Tank external dimensions: 535x335x265mm WxDxH
- 3-year warranty

**Special price - Ask us for a stocked item!**

**Special  
Price!**



Description	Code	Price (ex GST)	Sale Price (ex GST)
ST26 Stainless Steel Tank, 26L and T100 Immersion Thermostat, Temperature Range: Ambient +5°C to 100°C (or 0°C to 100°C with Cooling Accessory)	GRANT100-ST26	\$2,686.00	<b>\$2,549.00</b>
Lid Gabled Stainless Steel with Hinges for ST18 and ST26 Baths	GRANSTL26	\$597.00	<b>\$567.00</b>

**Bio-Strategy Ltd.**

T: 0800 34 24 66  
sales.nz@bio-strategy.com  
www.bio-strategy.com  
shop.bio-strategy.com