

Gyrozen Centrifuge on Sale!

VALID UNTIL:
31 December 2020

Ref. Code:
GYRO-CENT

Gyrozen offers a well-proven centrifuge design for Clinical Applications. Gyrozen guarantees the highest possible laboratory standards that you require for your daily work with industry-leading safety and quality.

- Imbalance cut-off
- Lid drop protection
- Automatic lid opening
- Very quiet operation
- Safety Lid Lock
- 10-100 programmes memory
- Automatic RPM/RCF conversion
- Autoclavable and robust rotors
- 2-year warranty
- We have them in stock!



624R

What is your application?

How many samples do you usually process at the same time?

Great Savings



406



416



1580R

Model	406	416/624R	1580/1580R
Description	The compact low-speed model, 406 is ideal for routine work with small sample sizes.	The 416 and 624R (refrigerated model) are compact bench top centrifuges. Big enough to accommodate 32 x vacutainers or 12 x 15ml tubes with extremely small footprint.	The 1580 and 1580R multi-purpose centrifuges are designed for large-volume daily working laboratories to save time with streamlining processes.
Capacity:	6 x 15ml tubes per run	32 x vacutainers / 12 x 15ml tubes / 4 x 50ml tubes	76 x 15ml tubes / 28 x 50ml tubes
Temperature Range:	N/A	624R (-10 to +40°C); 416 – N/A	1580R (-20 to +40°C); 1580 – N/A
Max. Speed:	4,000rpm / 2,075xg	4,800rpm / 3,888xg (624R) and 4,000rpm / 2,826xg (416)	15,000rpm / 24,249xg for fixed angle rotor and 5,000rpm / 5,394xg (swing out).
Applications:	Cell Separation, Blood Fractionation; Separation Work of Blood, Urine, And Stool; Sperm Washing.	General Laboratory Applications including; Biological Sample Separations of Cellular Materials, Blood, Urine, Sperm, etc.	Ideal for Cell Culture Applications, Clinical Preparation, Molecular Biology Work etc.

Description	Code	Price (ex GST)	Sale Price (ex GST)
406 Low-Speed General Centrifuge. Package includes: Centrifuge, Angle Rotor GRA-15-6A with Sleeves to Process 6 x 15ml Tubes. Adaptors for 3,4,9ml Tubes available as an option.	GYROGZ-0406	\$2,388.00	\$1,967.00
416 Low-Speed General Centrifuge. Package includes: Centrifuge, Four-Wing Swing-Out Rotor (GRS-S-100-4), 4 x Biosealed 100ml Buckets, 1 set of Adaptors of your choice to spin either 12 x 15ml Conical Tubes or 4 x 50ml Conical Tubes. Adaptors for 1.5 -100ml Tubes available as an option.	GYROGZ-0416	\$6,352.00	\$4,995.00
Low-Speed Refrigerated Centrifuge, Black. Package includes: Centrifuge, Four-Wing Swing-Out Rotor (GRS-SR-100-4), 4 x Biosealed 100ml Buckets, 1 set of Adaptors of your choice to spin either 12 x 15ml Conical Tubes or 4 x 50ml Conical Tubes. Adaptors for 1.5 -100ml Tubes available as an option.	GYROGZ-0624R-B	\$10,872.00	\$8,379.00
1580 High-Speed Multi-Purpose Centrifuge. Package includes: Centrifuge, Four-Wing Swing-Out Rotor (GRS-G-750-4), 4 x Biosealed 750ml Buckets, 1 set of Adaptors of your choice to spin either 56 x 15ml Conical Tubes or 20 x 50ml Conical Tubes. Adaptors for 3 -750ml Tubes available as an option.	GYROGZ-1580	\$13,612.00	\$10,995.00
1580R High-Speed Refrigerated Multi-Purpose Centrifuge. Package includes: Centrifuge, Four-Wing Swing-Out Rotor (GRS-G-750-4), 4 x Biosealed 750ml Buckets, 1 set of Adaptors of your choice to spin either 56 x 15ml Conical Tubes or 20 x 50ml Conical Tubes. Adaptors for 3 -750ml Tubes available as an option.	GYROGZ-1580R	\$17,823.00	\$14,675.00

Bio-Strategy Ltd.

T: 0800 34 24 66
sales.nz@bio-strategy.com
www.bio-strategy.com
shop.bio-strategy.com