

# Heidolph Vortex Mixers

## Reax Top / Reax Control Vortex Mixer

- The standard or precision model
- The shaking orbit of 5mm reliably and quickly achieves an even distribution
- Reax control with electronic speed control – the speed remains constant even in the low range and at load changes

## Multi Reax

- Shaking up to 26 samples simultaneously and achieving excellent mixing results
- Specific adaptors are not required
- Excellent mixing results are achieved with the 3mm shaking orbit, even with large samples with solid content
- Speed adjustment from 150 to 2,000rpm
- Timer function up to 999 minutes

## Incredibly strong and robust, the strong shaking motion guarantees excellent mixing results without exception.

- Designed and made in Germany: 3-year warranty on all devices and an average operational lifespan of 10-years
- A worthwhile investment: maintenance and spark-free motors and are excellently suited for years of continuous use
- The sealed housing reliably protects against corrosion
- Fast, even distribution, even with solids and highly viscous media: ideally suited for short- or long-term operation



VALID UNTIL:  
31 December 2020

Ref. Code:  
HEID-MIX



**In Stock!**  
**3-Year Warranty!**  
**Lifespan of 10-Years!**



**MADE IN GERMANY**

3-year warranty on all devices and an average operational lifespan of 10 years

Multi-stage quality checks in development and production

Premium service according to the "research made easy" principle

Description	Code	Price (ex GST)	Sale Price (ex GST)
Reax Top with Small Pick-Up Attachment	HEID541-10000-00	\$585.25	\$385.00
Reax Control with Small Pick-Up Attachment	HEID541-11000-00	\$861.12	\$689.00
Multi Reax	HEID545-10000-00	\$2,714.20	\$2,388.00

**Bio-Strategy Ltd.**

T: 0800 34 24 66  
sales.nz@bio-strategy.com  
www.bio-strategy.com  
shop.bio-strategy.com