

Receive a 25% discount on a new SpectraMax ABS or ABS Plus Reader

VALID UNTIL:
21 December 2020

Ref. Code:
MDEV-ABS

Affordable monochromator-based plate readers.

Run various assays within the Vis/UV-Vis range such as ELISA, Bradford, BCA, microbial growth, protein quantitation and more with the new small but mighty ABS and ABS Plus Absorbance readers monochromator-based plate readers.

- ABS reader works with 96- and 384-well plates in the VIS range (340-850nm)
- ABS Plus works with 96-, 384 well plates as well as cuvettes in the UV/VIS range (190-1000nm)



BENEFITS of the ABS and ABS PLUS Readers:

- Flexibility to choose your wavelength, unlike filter-based readers
- Elimination of well-to-well pipetting differences with Patented PathCheck[®] Sensor
- Preconfigured protocols such as ELISA and Bradford in SoftMax[®] Pro Software
- Easy and fast assay setup and data analysis

25% OFF, Order Today!



Bio-Strategy Ltd.

T: 0800 34 24 66
sales.nz@bio-strategy.com
www.bio-strategy.com
shop.bio-strategy.com