

Miele PG8583 Laboratory Washer Package

VALID UNTIL:
31 December 2020

Ref. Code:
MIEL-PG8583

\$1,000.00 discount, free detergent and neutraliser (value \$248.00), free shipping (Metro only) and commissioning (value >\$800.00).

Miele Professional PG8583 have been the standard in Laboratory Underbench/Free-Standing Washers for the past four years. Their stainless steel interior combining powerful 500l/min pump with integrated heating elements enable safe washing cycles up to 93°C. The EcoDry passive drying function automatically opens the door at the end of the run when temperature has dropped below 70°C allowing the load to dry and cool down faster. The unit also features spray arm and spray pressure monitoring with alerts to advise you if a blockage or pressure deviation has occurred.

- 15 programs including new "Mini" program for lightly soiled glassware
- ADP model including pump for AD water and waterproofing kit
- Dispensing pumps for liquid detergent and liquid neutraliser
- Compact: 600mm x 600mm x 850mm (820mm with lid removed) for under-bench use
- Many other options and accessories to select from

Package 1 includes upper and lower baskets and inserts/modules for washing wide neck and narrow neck glassware.



Package Deal!

Description	Code	Price (ex GST)	Sale Price (ex GST)
PG8583 Washer + Upper Basket and Inserts, Lower Basket for Modules, 2 Injector Modules and Waterproofing Kit.	PG8583 Package 1.	\$18,094.00	\$16,285.00

Please note: Price includes washer, waterproofing kit, A102 Upper Basket with 2 x AK12 inserts for wide neck glassware, A150 Lower Basket for modules and 2 x A301/5 modules (each with 3 x 6 = 18 injector nozzles; 160mm high), commissioning and training (customer to organise local plumber and electrician for installation and cover local freight).

Bio-Strategy Ltd.

T: 0800 34 24 66
sales.nz@bio-strategy.com
www.bio-strategy.com
shop.bio-strategy.com