

PHCbi MCO-170AICL-PE CO₂ Incubator

VALID UNTIL:
31 December 2020

Ref. Code:
SANY-CO2

The MCO-170AICL-PE CO₂ Incubator enhances your cell growth via an intelligent communication interface ensuring precise temperature and CO₂ control. The inCu-saFe copper enriched stainless-steel interior provides continuous germicidal protection against airborne contaminants on work surfaces. The ceramic based single beam dual detector infrared system provides continuous auto zero CO₂ calibration. The tempered glass sealed inner door allows easy viewing of your cultures without affecting the chamber environment. Downloading data is simply a matter of plugging in a USB stick. Supplied with 4 adjustable shelves (10 max.).

- 165 litres; ambient + 5 to 50°C; 0% to 20% CO₂
- Intelligent - easy to read colour touch-screen panel
- Data download via USB port
- Pro-active inCu-saFe decontamination
- Integrated side wall tray supports - less parts to remove for cleaning
- DHA air jacket provides stable temperature control and single beam dual detector infrared system provides precise CO₂ control



10% OFF!



| Description | Code | Price (ex GST) | Sale Price (ex GST) |
|---|---------------------|----------------|---------------------|
| CO ₂ Incubator, 165L, inCu-saFE, IR Dual Sensor and 4 Adjustable Shelves | SANYMCO-170CAICL-PE | \$9,970.00 | \$8,975.00 |

Bio-Strategy Ltd.

T: 0800 34 24 66
sales.nz@bio-strategy.com
www.bio-strategy.com
shop.bio-strategy.com