

PHCbi MPR-S313-PE Pharmaceutical Refrigerator

VALID UNTIL:
31 December 2020

Ref. Code:
SANY-PHARM

3-year warranty and including free remote alarm terminal with battery back up (MPR-48B1-PW valued at \$325.00).

The PHCbi MPR-S313-PE Pharmaceutical Refrigerator provides reliable storage for pharmaceuticals, samples and reagents. The stainless-steel interior is easy to clean, and the 5 adjustable clear coated hard steel wire shelves allow you to optimise your storage requirements. Featuring two defrost systems: cycle-defrost with no impact on chamber temperature and defrosting by defrost sensor where, if cycle defrost is not enough to remove ice due to frequent door openings, high ambient temperatures or moist loads, defrosting will occur automatically. The fan forced air circulation ensures precise temperature control and uniformity throughout the chamber.

- 2 to 14°C
- 340L capacity with compact footprint (ideal for narrow corridors)
- Double glazing glass sliding doors
- 15% energy saving with new cycle defrosting system that also prevents temperature spikes
- LED interior lighting linked to door opening/closing for on/off operation
- Door open alarm and high/low temperature alarms

**20%
OFF!**



Description	Code	Price (ex GST)	Sale Price (ex GST)
Pharmaceutical Refrigerator, 340L with Compact Footprint and Sliding Doors	SANYMPR-S313-PE	\$6,146.00	\$4,917.00

Bio-Strategy Ltd.

T: 0800 34 24 66
sales.nz@bio-strategy.com
www.bio-strategy.com
shop.bio-strategy.com